Sarstedt TC-Inserts



Einleitung

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Introduction

Sarstedt offers a broad range of cell culture inserts for 6, 12 and 24 well TC plates. The Inserts are made of a transparent polystyrene housing to which a polyethylene terephthalate (PET) track-etched membrane has been sealed. Both sides of the membrane are tissue culture-treated (TC) for improved cell adhesion. The special geometry of the hanging Sarstedt TC inserts suspends the insert in an off-center position within the culture plate well, facilitating convenient pipetting with e.g. serological pipettes by creating a larger space on one side (see Fig. 1 ①).

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Additionally, the lowered top edge of the PS housing facilitates optimal gas exchange (see Fig. 1 ②).

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Certification

- Non-pyrogenic/endotoxin-free
- Non-cytotoxic
- Individually, sterile wrapped in a blister pack

General information on membrane application fields:

- Transparent membranes are suitable for both light and electron microscopy.
- Translucent membranes are recommended for electron microscopy and TEER (transepithelial electrical resistance) experiments.
- Translucent membranes with a pore diameter of 0.4 μm allow, due to the high pore density, ideal basolateral diffusion needed for transport, secretion, diffusion and cytotoxicity studies.
- Membranes with small pore sizes are suitable for co-culture experiments.
- Membranes with large pore sizes are recommended for chemotactic, invasion and migration experiments.

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Instructions for use

- Peel off cover sheet from blister package containing the Sarstedt TC insert under aseptic conditions. Use sterile forceps to remove insert from the blister and place it into the well of a suitable TC plate. Load as many wells as needed for experiments and controls.
- Add pre-warmed culture medium to the well. Refer to table “Product and ordering information” for recommended working volumes.
- It is recommended to allow the membrane of each insert to become moistened with the culture medium by incubating them for several minutes.
- Add cells and culture medium to the insert and follow standard tissue culture incubation and feeding procedures.
- To remove cells from the membrane, follow your standard trypsinization or scraping protocol.
- For fixation and staining of cells, a broad spectrum of standard techniques can be performed. The membrane can be removed from the housing using a scalpel. The detached membrane does not curl up, simplifying further processing.

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**For further information, please visit
www.sarstedt.com**

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SARSTEDT AG & Co.