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

Disposable Vacuum Filter System

The design of SARSTEDT's disposable vacuum filter system features convenient handling for the sterilization or clarification of aqueous solutions in the laboratory. The PES (polyethersulfone) membrane is ideal for filtering culture media, biological fluids and other media components due to its short filtration time and low protein binding characteristics.

Filter units are available in 250 ml, 500 ml, and 1 L sizes with 0.22 µm pore size. Further, the 500 ml and 1 L filter units are available with 0.45 µm pore size. The 500 ml size is also available with 0.1 µm pore size.

Instruction For Use

Remove filter unit from bag. Check to insure that the lower reservoir has not loosened during shipping. Tighten if necessary.

- 1 + 2** Attach vacuum hose to the connector **before** coupling the connector to the vacuum unit. Ensure that the connector is plugged in tightly. Place unit on a flat surface such as a countertop and provide support if necessary. Remove the filter cover.
- 3** Pour solution into the filter funnel.
- 4** Replace the cover. Apply vacuum until solution is filtered or lower reservoir is full (at 50 kpa). Turn off vacuum and disconnect the hose.
- 5 + 6** Disconnect and discard the filter. If the lower reservoir is to be used as a storage bottle, the sterile cap supplied with the filter unit should be used to cap the bottle.

NOTE:

- Disposable vacuum filter units are intended for use by personnel familiar with safe laboratory practices.
- Eye protection is recommended when filter units are used under vacuum to protect against possible implosion injuries.
- Vessels not specifically designed for vacuum applications must not be substituted for receiver vessels due to the danger of implosion.
- Failure of the filter unit may result from improper use of vacuum.