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

# Disposable Bottle Top Filter

The design of SARSTEDT's disposable bottle top filter system features convenient handling for the sterilization or clarification of aqueous solutions in the laboratory. The PES (polyethersulfone) membrane is ideal for filtering culture media, biological fluids and other media components due to its short filtration time and low protein binding characteristics.

Bottle Top Filters are available in 250ml and 500ml sizes with 0.22 µm pore size. Further, the 500ml Bottle Top Filter is available with 0.45 µm pore size.

## Instruction For Use

Remove bottle top filter from bag.

- 1** Screw filter on neck of a bottle to be used as a reservoir. Neck of filter is designed for bottles with a GL45 thread size (45 mm).
- 2** Ensure that the connector is plugged in tightly. Place unit on a flat surface such as a countertop and provide support if necessary. Remove filter cover.
- 3** Pour solution into the filter funnel.
- 4** Replace the cover. Apply vacuum until solution is filtered or lower reservoir is full (at 50 kpa). Turn off vacuum and disconnect the hose. Disconnect and discard the filter.

## NOTE:

- Disposable vacuum filter units are intended for use by personnel familiar with safe laboratory practices.
- Eye protection is recommended when filter units are used under vacuum to protect against possible implosion injuries.
- Vessels not specifically designed for vacuum applications must not be substituted for receiver vessels due to the danger of implosion.
- Failure of the bottle top filter unit may result from improper use of vacuum.