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SARSTEDT Microvette® CB 200

DE	Gebrauchshinweis – SARSTEDT Microvette® CB 200	2 - 5
EN	Instructions for use – SARSTEDT Microvette® CB 200	6 - 9
BG	Инструкции за употреба – SARSTEDT Microvette® CB 200	10 - 13
CS	Návod k obsluze – SARSTEDT Microvette® CB 200	14 - 17
DA	Brugsanvisning – SARSTEDT Microvette® CB 200	18 - 21
EL	Οδηγίες χρήσεως – SARSTEDT Microvette® CB 200	22 - 25
ES	Instrucciones de uso – SARSTEDT Microvette® CB 200	26 - 29
ET	Kasutusjuhend – SARSTEDT Microvette® CB 200	30 - 33
FR	Mode d'emploi – SARSTEDT Microvette® CB 200	34 - 37
HR	Uputa za upotrebu – SARSTEDT Microvette® CB 200	38 - 41
HU	Használati utasítás – SARSTEDT Microvette® CB 200	42 - 45
IT	Istruzioni d'uso – SARSTEDT Microvette® CB 200	46 - 49
KO	사용 설명서 – SARSTEDT Microvette® CB 200	50 - 53
LT	Naudojimo instrukcijos – SARSTEDT „Microvette® CB 200“	54 - 57
LV	Lietošanas norādes – SARSTEDT Microvette® CB 200	58 - 61
NL	Gebruiksaanwijzing – SARSTEDT Microvette® CB 200	62 - 65
NO	Bruksanvisning – SARSTEDT Microvette® CB 200	66 - 69
PL	Instrukcja obsługi – próbówka Microvette® CB 200 SARSTEDT	70 - 73
PT	Instruções de utilização – SARSTEDT Microvette® CB 200	74 - 77
RO	Instrucțiuni de utilizare – SARSTEDT Microvette® CB 200	78 - 81
RU	Инструкция по применению – SARSTEDT Microvette® CB 200	82 - 85
SK	Návod na Použitie – SARSTEDT Microvette® CB 200	86 - 89
SL	Navodila za uporabo – Microvette® CB 200 SARSTEDT	90 - 93
SV	Bruksanvisning – SARSTEDT Microvette® CB 200	94 - 97
TH	คำแนะนำในการใช้งาน – ระบบการเจาะเลือดจากเส้นเลือดฝอย Microvette® CB 200 ของ SARSTEDT	98 - 101
TR	Kullanım için talimatlar – SARSTEDT Microvette® CB 200	102 - 105
ZH	使用说明 – SARSTEDT Microvette® CB 200	106 - 109

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Intended use

The Microvette® CB 200 is a blood collection system and is used with lancets to collect capillary blood using end-to-end capillaries. The Microvette® CB 200 is used for collecting, manual in-house transport and processing of capillary blood specimens in the clinical laboratory. The product is intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

The Microvette® CB 200 consists of a plastic receptacle, a 200 µl capillary with integrated holder and a colour-coded plastic push cap, and various additives (preparations). The selection of the additive depends on the analytical test method. It is specified by the manufacturer of the test reagents and/or the analysis equipment on which the test is conducted.

Microvette® screw cap colour codes:

Additive	Letter code	Cap colour based on ISO 6710	Cap colour based on BS 4851*
Serum tube			
Microvette® CB 200 with Serum	CAT	-	white
Microvette® CB 200 with Serum Gel	CAT	-	brown
EDTA tube			
Microvette® CB 200 with K3 EDTA**	K3E	violet	red
Tube with glycolysis inhibitors			
Microvette® CB 200 with Fluoride/K3 EDTA	FE	-	yellow

*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry.

**also available with soft cap with loop

SARSTEDT Microvette® CB 200 Serum/Serum Gel

The Microvette® with serum contains a clotting activator and is used for collecting serum. The Microvette® with serum gel contains a clotting activator and is used for collecting serum. After centrifugation, the gel forms a layer separating the coagulum and the serum. Microvettes with serum / serum gel are intended for clinical chemistry tests from serum.

SARSTEDT Microvette® CB 200 K3 EDTA

The Microvette® with K3 EDTA is coated with K3 EDTA as an anticoagulant and is intended for haematological tests.

SARSTEDT Microvette® CB 200 Fluoride/K3 EDTA

The Microvette® with fluoride/K3 EDTA is coated with K3 EDTA as an anticoagulant and contains fluoride for glycolysis inhibition.

NOTE: Fluoride (used for measuring glucose) can cause an increase in haemolysis. Further information on substances that can cause interference can be found in the respective instructions for use of the assay manufacturer.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils according to the guidelines and procedures in your facility. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects for capillary blood collection in appropriate disposal containers.
4. Under- or overfilling the Microvettes leads to an incorrect ratio of blood to the preparation/additive and can result in false analysis results.
5. The blood collected and processed with the Microvette® is not intended to be re-injected into the human body.
6. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year as indicated.

Storage

Store the product at room temperature.

Limitations

1. If blood samples are stored in the Microvette®, the stability of the analytes should be evaluated by the respective laboratory or taken from the instructions for use of the analysis device manufacturer.
2. Even if plasma or serum is separated by centrifuging the Microvettes and/or a barrier is present, all cells are not necessarily completely separated. Residual metabolism or natural degradation can affect the concentrations of analytes. The stability of the analyte should be evaluated with respect to the storage container and the conditions of the respective laboratory.
3. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use of the assay/analysis device.

Specimen collection and handling

READ ALL OF THIS DOCUMENT BEFORE STARTING TO DRAW CAPILLARY BLOOD.

Preparation for drawing capillary blood and materials required:

1. All necessary Microvettes, labelled according to size and additive.
2. Gloves, gown, goggles or other suitable protective gear to protect against pathogens or potentially infectious materials transmitted by blood.
3. Labels for specimen identification.
4. (Safety) lancet or incision lancet.
5. Disinfectant for cleaning the puncture site (follow the facility's guidelines for preparing the puncture site for taking blood samples). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

Recommended order of collection (CLSI* GP42):

1. Tubes with EDTA
2. Tubes with lithium heparin with/without gel
3. Tubes with glycolysis inhibitors
4. Tubes for coagulation tests of serum with/without gel

NOTE: The first drop of blood is discarded.


NOTE: Follow the regulations of your facility for the order of collection.

Capillary technique with the end-to-end capillary

General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING CAPILLARY BLOOD COLLECTION.

1. Select a suitable Microvette® for the required specimen.
2. Select a site for the puncture with a (safety) lancet or incision lancet. Please follow the guidelines of the facility for preparing the puncture site.

1. Open the Microvette® and save the stopper for later use. Put on the enclosed cap with the attached end-to-end capillary.
2. Hold the Microvette® horizontally or tilted slightly and collect the drops of blood with the end-to-end capillary.
3. When the capillary is entirely filled, allow the blood to flow vertically into the receptacle until the capillary tube is completely empty. Discard the cap with the end-to-end capillary and close the receptacle by pressing the stopper back into it.
4. Invert the specimen carefully.

Information on coagulation/on the Microvette® CB 200 Serum/Serum Gel:

After collecting blood, let the Microvette® Serum/Serum Gel coagulate for 30 minutes in an upright position to ensure a clear separation layer after centrifugation.

The recommended times are based on an intact coagulation process. Blood from patients with coagulation disorders due to their illness or under anticoagulant therapy requires more time for complete coagulation.

Centrifugation

WARNING!
Centrifuging Microvettes with cracks or centrifuging at an excessive centrifugal force can lead to breakage of the Microvettes, which can release potentially hazardous substances.

Please see our centrifugation recommendations for the various preparations in the table below.
Centrifuge inserts are to be selected according to the size of the Microvettes used. The relative centrifugal force has the following ratio to the selected rpm:

$RCF (g \text{ force}) = 11.2 \times r \times (\text{rpm}/1000)^2$,
"RCF": "Relative centrifugal force"
"rpm": "Revolutions per minute"
"r" [in cm]: "Radius of rotation" from the middle of the centrifuge to the bottom of the Microvette®.

Microvettes without gel can be centrifuged in centrifuges with a fixed-angle or swing-out rotor.





Microvettes with gel are designed solely for use in centrifuges with swing-out rotor. Centrifugation in centrifuges with a fixed-angle rotor is not validated by SARSTEDT and is not recommended. In addition, centrifugation should be conducted no later than 2 hours after collection, as this could otherwise yield false results.

Microvettes should be centrifuged in accordance with the centrifugation conditions listed below. If other conditions are used, they must be validated by the user.

It must be ensured that the Microvettes fit well in the centrifuge inserts. Microvettes that protrude beyond the insert can catch on the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge instructions for use for further information.

CAUTION! Do not remove broken Microvettes by hand.

Refer to the centrifuge instructions for use for information on disinfecting the centrifuge.

Image	Additive	Times and g values (=RCF*)				
		Standard recommendation		Alternative recommendation		
	Serum	5 min	10,000 x g	or	10 min	2,000 - 10,000 x g
	EDTA ²		2,000 x g			
	Glucose fluoride		10,000 x g			
	Serum gel ¹		4,000 - 10,000 x g			

¹ For tubes prepared with gel, we recommend using swing-out rotors.

² Only for plasma separation, not for haematological tests.

*RCF: Relative centrifugal force

Centrifugation at 20 °C

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled capillary blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Product-specific standards and guidelines as amended/updated

CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

WHO/DIL/LAB/99.1 Rev.02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

*CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Store away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
 CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.
 CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
 WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".
 *CLSI (Clinical and Laboratory Standards Institute)

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CLSI* GP39 'Tubes and Additives for Venous and Capillary Blood Specimen Collection', Approved Standard.

CLSI* GP42 'Collection of Capillary Blood Specimens', Approved Standard.

CLSI* GP44 'Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests', Approved Guideline.

WHO/DIL/LAB/99.1 Rev.02 'WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002'.

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CLSI* GP39 „Tubes and Additives for Venous and Capillary Blood Specimen Collection”, Approved Standard.

CLSI* GP42 „Collection of Capillary Blood Specimens”, Approved Standard.

CLSI* GP44 „Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests”, Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 „WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002”.

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CLSI* GP39 " Tubes and Additives for Venous and Capillary Blood Specimen Collection" , Approved Standard.

CLSI* GP42 " Collection of Capillary Blood Specimens" , Approved Standard.

CLSI* GP44 " Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests" , Approved Guideline.

WHO/DIL/LAB/99.1 Rev.02 " WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002" .

*CLSI (Clinical and Laboratory Standards Institute, 临床实验室标准化协会)

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