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SARSTEDT Microvette® CB 300

DE	Gebrauchshinweis – SARSTEDT Microvette® CB 300	2 - 5
EN	Instructions for use – SARSTEDT Microvette® CB 300	6 - 9
BG	Инструкции за употреба – SARSTEDT Microvette® CB 300	10 - 13
CS	Návod k obsluze – SARSTEDT Microvette® CB 300	14 - 17
DA	Brugsanvisning – SARSTEDT Microvette® CB 300	18 - 21
EL	Οδηγίες χρήσεως – SARSTEDT Microvette® CB 300	22 - 25
ES	Instrucciones de uso – SARSTEDT Microvette® CB 300	26 - 29
ET	Kasutusjuhend – SARSTEDT Microvette® CB 300	30 - 33
FR	Mode d'emploi – SARSTEDT Microvette® CB 300	34 - 37
HR	Uputa za upotrebu – SARSTEDT Microvette® CB 300	38 - 41
HU	Használati utasítás – SARSTEDT Microvette® CB 300	42 - 45
IT	Istruzioni d'uso – Microvette® CB 300 SARSTEDT	46 - 49
KO	사용 설명서 – SARSTEDT Microvette® CB 300	50 - 53
LT	Naudojimo instrukcijos – SARSTEDT „Microvette® CB 300“	54 - 57
LV	Lietošanas norādes – SARSTEDT Microvette® CB 300	58 - 61
NL	Gebruiksaanwijzing – SARSTEDT Microvette® CB 300	62 - 65
NO	Bruksanvisning – SARSTEDT Microvette® CB 300	66 - 69
PL	Instrukcja obsługi – próbówka Microvette® CB 300 SARSTEDT	70 - 73
PT	Instruções de utilização – Microvette® CB 300 SARSTEDT	74 - 77
RO	Instrucţiuni de utilizare – SARSTEDT Microvette® CB 300	78 - 81
RU	Инструкция по применению – SARSTEDT Microvette® CB 300	82 - 85
SK	Návod na Použitie – SARSTEDT Microvette® CB 300	86 - 89
SL	Navodila za uporabo – Microvette® CB 300 SARSTEDT	90 - 93
SV	Bruksanvisning – SARSTEDT Microvette® CB 300	94 - 97
TH	คำแนะนำในการใช้งาน – ระบบการเจาะเลือดจากเส้นเลือดฝอย Microvette® CB 300 ของ SARSTEDT	98 - 101
TR	Kullanım için talimatlar – SARSTEDT Microvette® CB 300	102 - 105
ZH	使用说明 – SARSTEDT Microvette® CB 300	106 - 109

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Intended use

The Microvette® CB 300 is a blood collection system and is used with lancets to collect capillary blood using End-to-End capillaries. The Microvette® CB 300 is used for collecting, manual in-house transport and processing of capillary blood specimens in the clinical laboratory. The product is intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

The Microvette® CB 300 consists of a plastic receptacle, a small stopper for closing the capillary opening and a colour-coded soft cap with loop, and various additives (preparations). The selection of the additive depends on the analytical test method. The method is specified by the manufacturer of the test reagents and/or the analysis equipment on which the test is conducted.

Microvette® CB 300 soft cap colour codes:

Additive	Letter code	Cap colour based on ISO 6710	Cap colour based on BS 4851*
Serum tube			
Microvette® CB 300 with serum in a conical inner receptacle	CAT	red	white
Lithium heparin tube			
Microvette® CB 300 with lithium heparin in a conical inner receptacle**	LH	green	orange
EDTA tube			
Microvette® CB 300 with K2 EDTA in a cylindrical inner receptacle	K2E	violet	red
Tube with glycolysis inhibitors			
Microvette® CB 300 with fluoride/heparin in a conical inner receptacle	FH	-	yellow

*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry.
 **also available with brown outer container with light protection properties

SARSTEDT Microvette® CB 300 Serum

The Microvette® CB 300 with serum contains a clotting activator and is used for collecting serum. Microvettes CB 300 with serum are intended for clinical chemistry tests from serum.

SARSTEDT Microvette® CB 300 Lithium heparin

The Microvette® CB 300 with lithium heparin is coated with lithium heparin as an anticoagulant and is used for collecting plasma. Microvettes CB 300 with lithium heparin are intended for clinical chemistry tests from plasma.

SARSTEDT Microvette® CB 300 K2 EDTA

The Microvette® CB 300 with K2 EDTA is coated with the anticoagulant K2 EDTA and is intended for haematological tests.

SARSTEDT Microvette® CB 300 Fluoride/heparin

The Microvette® CB 300 with fluoride/heparin is coated with the anticoagulant heparin and contains fluoride for glycolysis inhibition. It is used to test glucose levels.

NOTE: Fluoride (used for measuring glucose) can cause an increase in haemolysis. Further information on substances that can cause interference can be found in the respective instructions for use of the assay manufacturer.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself from blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils according to the guidelines and procedures in your facility. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects for capillary blood collection in appropriate disposal containers.
4. Under- or overfilling the Microvettes leads to an incorrect ratio of blood to the preparation/additive and can result in false analysis results.
5. The blood collected and processed with the Microvette® is not intended to be re-injected into the human body.
6. Do not use the products after the expiry date. The expiry date is the last day of the month and year as indicated.

Storage

Store the product at room temperature.

Limitations

1. If blood samples are stored in the Microvette®, the stability of the analytes should be evaluated by the respective laboratory or taken from the instructions for use of the analysis device manufacturer.
2. Even if plasma or serum is separated by centrifuging the Microvettes and/or a barrier is present, all cells are not necessarily completely separated. Residual metabolism or natural degradation can affect the concentrations of analytes. The stability of the analyte should be evaluated with respect to the storage container and the conditions of the respective laboratory.
3. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use of the assay/analysis device.

Specimen collection and handling

READ ALL OF THIS DOCUMENT BEFORE STARTING TO DRAW CAPILLARY BLOOD.

Preparation for drawing capillary blood and material needed:

1. All necessary Microvettes, labelled by size and additive.
2. Gloves, gown, goggles or other suitable protective gear to protect against bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. (Safety) lancet or incision lancet.
5. Disinfectant for cleaning the puncture site (follow the facility's guidelines for preparing the puncture site for taking blood samples). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

Recommended order of collection (CLSI* GP42):

1. Tubes with EDTA
2. Tubes with lithium heparin with/without gel
3. Tubes with glycolysis inhibitors
4. Tubes for coagulation tests of serum with/without gel

NOTE: The first drop of blood is discarded.


NOTE: Follow the regulations of your facility for the order of collection.

Capillary technique with the End-to-End capillary

General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING CAPILLARY BLOOD COLLECTION.

1. Select a suitable Microvette® for the required specimen.
2. Select a site for the puncture with a (safety) lancet or incision lancet. Please follow the guidelines of the facility for preparing the puncture site.

- 1 Take the trumpet-shaped Microvette® tube out of the carrier tube.
- 2 Hold the Microvette® horizontally or tilted slightly and collect the drops of blood with the tip. The receptacle fills automatically due to capillary action.
- 3 After collecting the blood, affix the attached small stopper firmly onto the tip in horizontal or slightly tilted position.
- 4 Then push the top stopper on.
- 5 Replace the closed receptacle in the freestanding carrier tube.
- 6 Mix the specimen carefully and gently!

Information on coagulation/on the Microvette® CB 300 Serum:

After collecting blood, let the Microvette® CB 300 Serum coagulate for 30 minutes in an upright position to ensure a clear separation layer after centrifugation.

The recommended times are based on an intact coagulation process. Blood from patients with coagulation disorders due to their illness or under anticoagulant therapy requires more time for complete coagulation.

Centrifugation

WARNING!

Centrifuging Microvettes with cracks or centrifuging at an excessive centrifugal force can lead to breakage of the Microvettes, which can release potentially hazardous substances.

Please see our centrifugation recommendations for the various additives in the table below.

Centrifuge inserts are to be selected according to the size of the Microvettes used. The relative centrifugal force has the following ratio to the selected rpm:

$$RCF (g\text{-force}) = 11.2 \times r \times (\text{rpm}/1000)^2,$$

"RCF": "Relative centrifugal force"

"rpm": "Revolutions per minute" (min⁻¹),

"r" [in cm]: "Radius of rotation from the middle of the centrifuge to the bottom" of the Microvette®.

Microvettes without gel can be centrifuged in centrifuges with a fixed-angle or swing-out rotor.



Centrifugation in centrifuges with a fixed-angle rotor is not validated by SARSTEDT and is not recommended. In addition, centrifugation should be conducted no later than 2 hours after collection, as this could otherwise yield false results.

Microvettes should be centrifuged in accordance with the centrifugation conditions listed below. If other conditions are used, they must be validated by the user.

It must be ensured that the Microvettes fit well in the centrifuge inserts. Microvettes that protrude beyond the insert can catch on the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge instructions for use for further information.

CAUTION! Do not remove broken Microvettes by hand.

Refer to the centrifuge instructions for use for information on disinfecting the centrifuge.

Image	Additive	Times and g values (=RCF*)				
		Standard recommendation		or	Alternative recommendation	
	Serum	5 min	10,000 x g			10 min
	Lithium heparin		2,000 x g			
	Glucose fluoride					

*RCF: Relative centrifugal force

Centrifugation at 20°C

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled capillary blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Product-specific standards and guidelines as amended/updated

CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

WHO/DIL/LAB/99.1 Rev.02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

*CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

Article number

Batch number

Use by

CE marking

In vitro diagnostic device

Follow the instructions for use

If reused: Risk of contamination

Store away from sunlight

Store in a dry place

Manufacturer

Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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CLSI* GP39 „Tubes and Additives for Venous and Capillary Blood Specimen Collection”, Approved Standard.

CLSI* GP42 „Collection of Capillary Blood Specimens”, Approved Standard.

CLSI* GP44 „Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests”, Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 „WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002”.

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CLSI* GP39 „Tubes and Additives for Venous and Capillary Blood Specimen Collection“, Approved Standard.

CLSI* GP42 „Collection of Capillary Blood Specimens“, Approved Standard.

CLSI* GP44 „Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests“, Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 „WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002“.

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- CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
- CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.
- CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
- WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations Geneva: World Health Organization, 2002".
- *CLSI (Clinical and Laboratory Standards Institute)

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CLSI* GP39 " Tubes and Additives for Venous and Capillary Blood Specimen Collection" , Approved Standard.

CLSI* GP42 " Collection of Capillary Blood Specimens" , Approved Standard.

CLSI* GP44 " Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests" , Approved Guideline.

WHO/DIL/LAB/99.1 Rev.02 " WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002" .

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