S-Monovette® cfDNA Exact

DE	Gebrauchshinweis – S-Monovette® cfDNA Exact	3 - 6
EN	Instructions for use – S-Monovette® cfDNA Exact	7 - 10
BG	Инструкции за употреба – S-Monovette® cfDNA Exact	11 - 14
BR	Instruções de utilização – S-Monovette® cfDNA Exact	15 - 18
CS	Návod k obsluze – S-Monovette® cfDNA Exact	19 - 22
DA	Brugsanvisning – S-Monovette® cfDNA Exact	23 - 26
EL	Οδηγίες χρήσεως – S-Monovette® cfDNA Exact	27 - 30
ES	Instrucciones de uso – S-Monovette® cfDNA Exact	31 - 34
ET	Kasutusjuhend – S-Monovette® cfDNA Exact	35 - 38
FR	Mode d'emploi – S-Monovette® cfDNA Exact	39 - 42
HR	Uputa za upotrebu – S-Monovette® cfDNA Exact	43 - 46
HU	Használati utasítás – S-Monovette® cfDNA Exact	47 - 50
IT	Istruzioni d'uso – S-Monovette® cfDNA Exact	51 - 54
KO	사용 설명서 – S-Monovette® cfDNA Exact	55 - 58
LT	Naudojimo instrukcijos – S-Monovette® cfDNA Exact	59 - 62
LV	Lietošanas norādes – S-Monovette® cfDNA Exact	63 - 66
NL	Gebruiksaanwijzing – S-Monovette® cfDNA Exact	67 - 70
NO	Bruksanvisning – S-Monovette® cfDNA Exact	71 - 74
PL	Instrukcja obsługi – S-Monovette® cfDNA Exact	75 - 78
PT	Instruções de utilização – S-Monovette® cfDNA Exact	79 - 82
RO	Instrucțiuni de utilizare – S-Monovette® cfDNA Exact	83 - 86
RU	Инструкция по применению – S-Monovette® cfDNA Exact	87 - 90
SK	Návod na Použitie – S-Monovette® cfDNA Exact	91 - 94
SL	Navodila za uporabo – S-Monovette® cfDNA Exact	95 - 98
SV	Bruksanvisning – S-Monovette® cfDNA Exact	99 - 102
TH	คำแนะนำในการใช้งาน – S-Monovette® cfDNA Exact	103 - 106
TR	Kullanım için talimatlar – S-Monovette® cfDNA Exact	107 - 110
US	Instructions for use – S-Monovette® cfDNA Exact	111 - 114
ZH	使用说明 – S-Monovette® cfDNA Exact	115 - 118

S-Monovette® cfDNA Exact

DE

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Instructions for Use


Purpose

The S-Monovette® cfDNA Exact is used for venous blood collection with the aspiration or vacuum technique. The sample volume is 9.2 ml (of which 8.4 ml is blood).

The product is used to immediately stabilise the cfDNA present in the blood (preventing degradation of the cfDNA and the release of genomic DNA from nucleated cells). The S-Monovette® cfDNA Exact separates plasma from the stabilised blood using centrifugation. The S-Monovette® cfDNA Exact standardises preanalytics in a closed system while the sample is being stored and transported.

Required materials

To obtain venous blood and plasma for cfDNA analytics, you will need the following products, some of which are not included in the S-Monovette® cfDNA Exact scope of supply *:


- ① A tourniquet, e.g. REF 95.1006
- ② A Safety-Multify® needle (REF 85.1637.235, 85.1638.235, 85.1640.235 or 85.1642.235), length: 200 mm
- ③ An empty/discard tube (e.g. S-Monovette® Neutral or other S-Monovette®) to fill the tube of the Safety-Multify® needle
- ④ **S-Monovette® cfDNA Exact, REF 01.2040.001**
- ⑤ A centrifugation cap for plasma separation (can be ordered separately if required under REF 65.729)

Handling


The stabiliser solution pre-filled in the S-Monovette® cfDNA Exact is harmful to health. You must therefore follow the prescribed measures for blood collection and only use the materials prescribed in these instructions for use. Disinfect the puncture site beforehand as per the internal hygiene plan. Always fill the S-Monovette® cfDNA Exact as the last S-Monovette® when making multiple collections. Blood can be collected with the S-Monovette® cfDNA Exact using the aspiration or vacuum technique.

When using the aspiration technique, fill the S-Monovette® cfDNA Exact after connecting by slowly pulling up the plunger.


When using the vacuum technique, pull up and snap the plunger until it clicks into the rear snap-in position before connecting the S-Monovette® cfDNA Exact.


- 1) Put on tourniquet (do not keep in place for longer than 30 seconds). Puncture vein in a suitable place.


- 2) Connect empty/discard tube, fill Safety-Multify® needle with blood.


- 3) Loosen tourniquet, fill empty/discard tube and remove.

- 4) Hold the S-Monovette® cfDNA Exact vertically below arm level and connect. Fill the S-Monovette® cfDNA Exact by pulling up the plunger (aspiration technique).

- 5) Disconnect, pull the plunger into the snap-in position and snap (aspiration technique).

- 6) Mix preparation by inverting the S-Monovette® cfDNA Exact at least five times after collecting the blood.

Sending samples

Samples are sent as full blood in the blood collection tube. As samples are often collected in a different place, it is often necessary to transport the sample to the lab after blood collection. In order to ensure secure transport, SARSTEDT offers a complete packaging system for the S-Monovette®, which meets the requirements of packaging instruction P 650, the ADR, RID, ICAO and IATA:

- Protective container with absorbent lining (REF 78.898) and screw cap (REF 65.679) as secondary packaging
- Shipping box (REF 95.900, 95.901 & 95.902)
- Shipping packaging for refrigerated transport (REF 95.995, 95.1123, 95.1124 or 95.1625)
- Styrofoam casing for refrigerated transport (REF 95.1011)

Sample storage


In the S-Monovette® cfDNA Exact, the cfDNA in the blood sample remains stable for at least:

- 14 days at 4°C to 37°C

It is not possible to freeze the blood samples.

Plasma separation

After sample collection, the sample volume in the S-Monovette® cfDNA Exact is 9.2 ml (of which 8.4 ml is blood). Separating the plasma to obtain the cfDNA is a multi-stage centrifugation process, which can be carried out as follows for example:

1. Replace the red screw-on lid with the centrifugation cap included.
2. Centrifuge for 10 minutes at 2,000–3,000 x g at room temperature.
3. Transfer the plasma into reaction vessels (e.g. REF 72.706.200 (1.5 ml), 72.695.200 (2.0 ml) or 72.701.400 (5.0 ml)).
4. Centrifuge the plasma for 15 minutes at 15,000 x g at room temperature.
5. Transfer the plasma into new vessels to isolate the cfDNA and/or to store (-80°C) until isolation.

Suitable products for cfDNA isolation are:

1. Manual purification systems

- NucleoSpin® Plasma XS-Kit, Macherey-Nagel, REF 740900.50
- NucleoSnap® DNA Plasma-Kit, Macherey-Nagel, REF 740300.50
- QIAamp® Circulating Nucleic Acid Kit, Qiagen, REF 55114
- MagMAX™ Cell-Free DNA Isolation Kit, ThermoFisher Scientific, REF A29319

2. Automated purification systems

- Chemagic™ cfDNA 5k Kit H24, Perkin Elmer, REF CMG-1304
- InviMag® Free Circulating DNA Kit/IG, Invitex Molecular, REF 2439320400

You can find more information on the SARSTEDT website.

Warnings

- Do not loosen the tourniquet while the S-Monovette® cfDNA Exact is connected.
- Blood and blood components that are collected and processed in the S-Monovette® cfDNA Exact are not suitable for infusion or introduction into the human body.
- The usual precautionary measures should be observed when handling laboratory reagents.
- Only use S-Monovettes for blood collection and not for injection.

Warning

Hazard statements

H317 May cause an allergic skin reaction.

Precautionary statements

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
 P272 Contaminated work clothing should not be allowed out of the workplace.
 P280 Wear protective gloves/protective clothing/eye protection/face protection.
 P302+P352 IF ON SKIN: Wash with plenty of water.
 P321 Specific treatment (see remarks on this label).
 P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
 P362+P364 Take off contaminated clothing and wash it before reuse.
 P501 Dispose of contents/container to waste disposal site according to legal regulations.

Safety data sheet at www.sarstedt.com

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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


- 2) Connect empty/discard tube, fill Safety-Multitype® needle with blood.

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Warning

Hazard statements

May cause an allergic skin reaction

Precautionary statements

Avoid breathing dust/fume/gas/mist/vapors/spray.
 Contaminated work clothing must not be allowed out of the workplace.
 Wear protective gloves/protective clothing/eye protection/face protection.
 If on skin: Wash with plenty of water.
 Specific treatment (see remarks on this label).
 If skin irritation or rash occurs: Get medical advice/attention.
 Wash contaminated clothing before reuse.
 Dispose of contents/container to waste disposal site according to legal regulations.

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SARSTEDT AG & Co. KG