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SARSTEDT Zwischenbodenröhre CSF

DE	Gebrauchshinweis – SARSTEDT Zwischenbodenröhre CSF	2 - 4
EN	Instructions for use – SARSTEDT CSF false bottom tube	5 - 7
BG	Инструкции за употреба – SARSTEDT Епруветка с междинно дъно ЦСТ	8 - 10
CS	Návod k obsluze – Zkumavka s falešným dnem CSF SARSTEDT	11 - 13
DA	Brugsanvisning – SARSTEDT Mellembundsrør CSF	14 - 16
EL	Οδηγίες χρήσεως – Σωληνάριο ENY με ενδιάμεσο πάτο SARSTEDT	17 - 19
ES	Instrucciones de uso – Tubo con fondo intermedio CSF SARSTEDT	20 - 22
ET	Kasutusjuhend – SARSTEDT Zwischenbodenröhre CSF	23 - 25
FR	Mode d'emploi – SARSTEDT Tubes à fond intermédiaire LCS	26 - 28
HR	Upute za uporabu – SARSTEDT epruveta CSF s lažnim dnom	29 - 31
HU	Használati utasítás – SARSTEDT köztes aljú CSF mintavételi csővek	32 - 34
IT	Istruzioni d'uso – SARSTEDT Provetta con doppio fondo LCR	35 - 37
KO	사용 설명서 – SARSTEDT CSF 이중바닥 튜브	38 - 40
LT	Naudojimo instrukcijos – SARSTEDT Zwischenbodenröhre CSF	41 - 43
LV	Lietošanas norādes – SARSTEDT papilddibena mēģene CSF	44 - 46
NL	Gebruiksaanwijzing – SARSTEDT CSF-buisje met tussenbodem	47 - 49
NO	Bruksanvisning – SARSTEDT false bottom rør CSF	50 - 52
PL	Instrukcja obsługi – SARSTEDT Probówka z podwójnym dnem na PMR	53 - 55
PT	Instruções de utilização – Tubo com fundo falso LCR SARSTEDT	56 - 58
RO	Instrucțiuni de utilizare – Tub cu fund conic pentru LCR SARSTEDT	59 - 61
RU	Инструкция по применению – Пробирка с двойным дном для спинномозговой жидкости (СМЖ) SARSTEDT	62 - 64
SK	Návod na Použitie – SARSTEDT Skúmavka CSF s falošným dnom	65 - 67
SL	Navodila za uporabo – epruveta CST z vmesnim dnom SARSTEDT	68 - 70
SV	Bruksanvisning – SARSTEDT Falsebottomrör för CSF	71 - 73
TH	คำแนะนำในการใช้งาน – หลอดเก็บตัวอย่าง CSF แบบมีพื้นคั่นของ SARSTEDT	74 - 76
TR	Kullanım için talimatlar – SARSTEDT BOS Ara Tabanlı Tüp	77 - 79
ZH	使用说明 – SARSTEDT带隔层样品管CSF	80 - 82

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
Systembezogene Normen:

DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes
 MDSAP: Medical Device Single Audit Program gemäß ISO 13485
 DIN EN ISO 14971: Medical devices - Application of risk management to medical devices
 IVDR (EU) 2017/746: Verordnung über *in-vitro*-Diagnostika

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Intended use

The CSF false bottom tube is a specimen receptacle for CSF (cerebrospinal fluid) which has been designed for the collection, transportation, processing and storage of CSF for *in vitro* diagnostics. The product is intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

The CSF false bottom tube consists of a clear plastic receptacle with a round bottom and an elevated conical base. The tube has a transparent plastic label. The nominal volume is 2.5 ml and is printed on the tube as a fill level with the additional information "Fill to mark". The product is individually sterile-packed with an enclosed blue screw cap.

Safety information and warnings

Do not use the product if the blister packaging is damaged.

1. General precautions: Wear gloves and other general personal protective equipment to prevent possible exposure to potentially infectious specimen material and transmitted pathogens.
2. Handle all biological specimens and sharp collection utensils according to the established guidelines and procedures in your facility. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Check each individual package for damage before use and do not use if packaging is damaged.
4. Dispose of all sharp objects (cannulae) for specimen collection in the appropriate hazardous biological waste containers.
5. The CSF false bottom tube and cap are intended for single use only. Dispose of all products in the hazardous biological waste containers.
6. Avoid under- or over-filling and foam formation.
7. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year as indicated.

Storage

Store the product at room temperature.

Transportation

The product is a primary receptacle according to the ADR (P650 packaging instruction) and IATA regulations.

Limitations

1. The storage period and temperature of a filled false bottom tube are dependent on the expiry date of the parameters for investigation.
2. Check the CSF specimen quality. Note the information provided by the test manufacturer regarding suitability for haemolysis and precipitates for clinical tests.

Collecting and handling specimens

READ THIS DOCUMENT IN FULL BEFORE COMMENCING THE PUNCTURE.

Equipment needed to collect a specimen:

1. Appropriate protective clothing and general equipment for CSF collection in line with the established hygiene guidelines at your facility.
2. All required specimen receptacles as per the request.
3. Labels for specimen identification.
4. Sharps disposal container for the safe disposal of used equipment.

N.B.: Check whether specific pre-analytical handling processes are necessary according to the test manufacturer and follow these steps.

COLLECT THE SPECIMEN IN LINE WITH THE ESTABLISHED GUIDELINES AND PROCEDURES IN YOUR FACILITY.

Handling the CSF false bottom tube

1. To use the sterile CSF false bottom tube, open the blister pack using the peel-off technique over a sterile cloth. The transparent film must be facing upwards.
N.B.: The blister pack contains a CSF false bottom tube and a screw cap which is not attached.
2. Allow the cerebrospinal fluid to drip directly into the CSF false bottom tube until it reaches the fill level mark.
3. Put the cap on and carefully close the CSF false bottom tube.
4. Place the CSF false bottom tube in an upright position according to the intended storage conditions for the specimen material.
5. Affix the specimen identification.

Centrifugation

WARNING!

Use suitable tube holders or inserts only. Centrifugation of a cracked CSF false bottom tube, or centrifugation using excessive centrifugal force can cause breakage of the CSF false bottom tube, releasing potentially hazardous substances.

Centrifuge adapters must be selected according to the size of the respective CSF false bottom tubes. The relative centrifugal force has the following ratio to the selected rpm:

$$\text{RCF (G-force)} = 11.2 \times r \times (\text{rpm}/1000)^2,$$

"RCF": "Relative centrifugal force"

"RPM" or "rpm": n = "Revolutions per minute"

"r" [in cm]: "Radius of rotation from the centre of the centrifuge to the base of the elevated section of the CSF false bottom tube."

CSF false bottom tubes should be centrifuged in accordance with the centrifugation conditions set out below. If other conditions are used, they must be validated by the respective user.

Ensure that the CSF false bottom tube fits snugly in the centrifuge adapter. Tubes that protrude beyond the adapter can catch on the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge operating manual for further information.

CAUTION! Do not remove broken CSF false bottom tubes by hand.

Refer to the centrifuge operating manual for instructions on disinfecting the centrifuge.

Centrifugation recommendation:

The recommended RCF force (g-force) for the centrifugation of CSF specimens with precipitates is 2,000 x g for 10 minutes.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled CSF false bottom tubes must be disposed of in suitable hazardous biological waste containers which can subsequently be autoclaved and burned.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

US standards / ISO standards (as amended/updated)

Product-specific standards

DIN EN 14254: *In-vitro* diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans
 DIN EN 14401 Rigid plastics containers – Methods to test the effectiveness of closures

Standards for labelling:

DIN EN 1041: Information supplied by the manufacturer of medical devices
 DIN EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
 Commission Regulation (EU) No 207/2012: Regulation on electronic instructions for use of medical devices

Standards for packaging:

DIN EN ISO 11607: Packaging for terminally sterilized medical devices
 DIN EN 868-5: Packaging for terminally sterilized medical devices – Part 5 Sealable pouches and reels of porous materials and plastic film construction – Requirements and test methods

Standards for sterilization:

DIN EN ISO 11135: Sterilization of health-care products – Ethylene oxide – Requirements for the development, validation and routine control of a sterilization process for medical devices
 DIN EN ISO 11737-1: Sterilization of health-care products – Microbiological methods – Part 1: Determination of a population of microorganisms on products
 DIN EN 556-1: Sterilization of medical devices – Requirements for medical devices to be designated "STERILE" – Part 1: Requirements for terminally sterilized medical devices


System-related standards:

DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes
 MDSAP: Medical Device Single Audit Program in accordance with ISO 13485
 DIN EN ISO 14971: Medical devices - Application of risk management to medical devices
 IVDR (EU) 2017/746: Regulation on *in vitro* diagnostic medical devices

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Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE mark
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	Risk of contamination if reused
	Keep away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

The following is also applicable to sterile products:

	Do not resterilise
	Ethylene oxide sterilisation
	Do not use if package is damaged
	Single sterile barrier system with external protective packaging

Technical modifications reserved

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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Systemrelaterede standarder:

DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes
MDSAP: Medical Device Single Audit Program til ISO 13485
DIN EN ISO 14971: Medical devices - Application of risk management to medical devices
IVDR (EU) 2017/746: Forordning om *in-vitro*-diagnostik

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Hansson O, Rutz S, et al. *Pre-analytical protocol formeasuring Alzheimer's disease biomarkers in fresh CSF*, DOI: 10.1002/dad2.12137

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Norme specifiche del prodotto:

DIN EN 14254: *In-vitro* diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans
DIN EN 14401 Rigid plastics containers - Methods to test the effectiveness of closures

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DIN EN 1041: Information supplied by the manufacturer of medical devices
DIN EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
(UE) 207/2012: Regolamento relativo alle istruzioni per l'uso elettroniche dei dispositivi medici

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DIN EN ISO 11607: Packaging for terminally sterilized medical devices
DIN EN 868-5: Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods

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DIN EN ISO 11135: Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices
DIN EN ISO 11737-1: Sterilization of health care products - Microbiological methods - Part 1: Determination of a population of microorganisms on products
DIN EN 556-1: Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices

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DIN EN ISO 13485: Medical devices - Quality management systems - Requirements for regulatory purposes
MDSAP: Medical Device Single Audit Program secondo ISO 13485
DIN EN ISO 14971: Medical devices - Application of risk management to medical devices
IVDR (UE) 2017/746: Regolamento relativo ai dispositivi medico-diagnostici *in vitro*

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Standarde privind marcajele:

DIN EN 1041: Information supplied by the manufacturer of medical devices
 DIN EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
 (EU) 207/2012: Regulament privind instrucțiunile electronice de utilizare a dispozitivelor medicale

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Standarde referitoare la sisteme:

DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes
 MDSAP: Medical Device Single Audit Program conform ISO 13485
 DIN EN ISO 14971: Medical devices - Application of risk management to medical devices
 IVDR (EU) 2017/746: Regulament privind diagnosticarea *in-vitro*

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Hygienemaßnahmen bei Liquorpunktionen, Liquorableitungen und Injektionen am ZNS. [Măsuri de igienă pentru punctțiile de lichid cefalorahidian, drenajul lichidului cefalorahidian și injecțiile la nivelul SNC.] Hyg Med 2011; 36 - 4, pag. 141-142
 Hansson O, Rutz S, et al. *Pre-analytical protocol formeasuring Alzheimer's disease biomarkers in fresh CSF*, DOI: 10.1002/dad2.12137
 Leitlinien für Diagnostik und Therapie in der Neurologie, Lumbalpunktion und Liquordiagnostik, Entwicklungsstufe: [Orientări pentru diagnosticare și terapie în neurologie, puncție lombară și diagnosticare LCR, etapa de realizare:] S1, număr de înregistrare AWMF: 030/141

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DIN EN 1041:Information supplied by the manufacturer of medical devices
DIN EN ISO 15223-1:Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1:General requirements (EU) 207/2012:医疗产品电子使用说明条例

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DIN EN ISO 13485:Medical devices - Quality management systems – Requirements for regulatory purposes
MDSAP:根据Medical Device Single Audit Program规定, 符合ISO 13485
DIN EN ISO 14971:Medical devices - Application of risk management to medical devices
IVDR (EU) 2017/746:关于体外诊断医疗设备的规定

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