

Differentiation of Mesenchymal Stem Cells within GrowDex[®] using SARSTEDT TC inserts

SARSTEDT AG & Co. KG



Authors

Jonathan Sheard^{1,2}, Ioannis Azoidis², Dariusz Widera²

¹Sheard BioTech Ltd, ²Stem Cell Biology and Regenerative Medicine Group, University of Reading, UK

Introduction

Mesenchymal stem cells (MSCs) are clinically relevant tripotent adult stem cells which give rise to bone (osteogenic), cartilage (chondrogenic) and fat (adipogenic) cells [1]. Since their discovery, differentiation of MSCs has been traditionally conducted on flat two-dimensional (2D) glass or polystyrene surfaces. 2D cell culture is known to result in unnatural cell polarity and morphology in addition to lacking the three-dimensional extracellular matrix. Moreover, MSCs are known to sensitively respond to topological and mechanical cues. Thus, differentiation in 3D can be considered more physiological than under 2D conditions.

Indeed, modulated osteogenic and adipogenic differentiation properties of MSCs have been reported [3, 4, 5]. In addition, delivery of bone forming cells within a known and controlled 3D scaffold may help towards regenerating bone fractures or critical size defects [6, 7].

Recently, we reported successful culture expansion of MSCs within three dimensional nanofibrillar cellulose (GrowDex[®]) [8], [application note 22].

Here, we describe osteogenic and adipogenic differentiation of human MSCs embedded within GrowDex[®] in SARSTEDT 24-well TC inserts.





Fig. 1
Images showing MSCs stained with Oil Red O following 21 days of incubation in either control (left) or adipogenic differentiation (right) media. Deposition of intracellular lipid droplets can be seen within differentiated cells (right). Scale bar is 50μm.

Materials

- Palatal adipose tissue-derived Mesenchymal Stem Cells obtained from adult donors with written informed consent with an approval from the Ethics Committee of Dental Faculty of Selcuk University (approval number 2012-08)
- Complete Media: DMEM (High Glucose) supplemented with 10 % FBS, 1 % Pen/Strep and 1 % L-Glutamine (Sigma Aldrich)
- Differentiation Media: StemPro™ Adipogenesis or Osteogenesis Differentiation Kit (Cat# A1007001 and Cat# A1007201, ThermoFisher) supplemented with 1 % Pen/Strep
- 24 well tissue culture inserts: 3.0 μm pore PET membrane (Cat# 83.3932.300, SARSTEDT)
- Low adhesion 24 well cell-culture plate (Cat# 83.3922.500, SARSTEDT)

- GrowDex®, 1.5 % (Cat# 100 103 005, UPM)
- 4 % Paraformaldehyde (PFA)
- Oil Red O: 300 mg Oil Red O in 100 ml 99% isopropanol (Cat# O0625-25G, Sigma Aldrich)
- Alizarin Red S: 2 g Alizarin Red S in 100 ml distilled water, pH 4.1 with 0.1 % NH₄OH, filtered (A5533-25G, Sigma Aldrich)
- EVOS Live imaging system (ThermoFisher)



Methods

1. Adipose derived and palatal derived MSCs were cultured in complete media and incubated in a humidified incubator at 37°C and 10 % CO₂.
2. Following trypsinization, cells were resuspended in complete media at a concentration of 1x10⁶ cells/ml.
3. Cells were mixed with the appropriate volume of complete media and GrowDex® to provide a final concentration of 0.2 % GrowDex® seeded with 2x10⁵ cells (2,000 cells/µl).
4. Embedded cells were transferred into 24-well tissue culture inserts and fed with 500 µl complete media within the outer well.
5. Medium was changed after 3 days to either complete media (control), or StemPro™ Adipogenic/Osteogenic media and incubated at 37°C with 10 % CO₂.

6. Control and differentiation media were replaced ever 2 to 3 days and cells were maintained for 21 days.
7. Following the differentiation period, media was removed from the wells and inserts and cells were washed with PBS and then fixed with 4% PFA for 30 mins.
8. For Oil Red O staining: 3 parts of the stock solution were mixed with 2 parts DI water and then filtered. Following fixation, cells were washed with sterile H₂O. Staining solution was added to the cells for 5 mins. For Alizarin red staining: following fixation, Alizarin red staining solution was added to cells and incubated for 45 mins at RT in the dark.
9. The staining solutions were removed and unbound dye was washed off by 3 to 5 washing steps with ddH₂O.
10. Images were taken using the EVOS imaging system.

Results

Following 21 days of differentiation treatments, MSCs showed lipid accumulation (Fig.1) and calcium deposition (Fig.2). Oil Red O staining of cells treated with adipogenic differentiation media revealed strong staining of intracellular lipids (Fig.1, right) in contrast to the control (Fig.1, left). Alizarin red staining of cells cultivated in

osteogenic differentiation media showed high levels of calcium deposition (Fig.2, right) compared to the control (Fig.2, left).

Conclusions

Adipogenic or osteogenic differentiation potential is an important identification criteria for MSCs [1, 2]. Compared to 2D, differentiation of MSCs in an 3D environment represents a more physiological approach and is known to enhance the levels of adipogenesis and osteogenesis [3, 4, 5]. Additionally, MSCs embedded within 3D scaffolds show a greater potential for osteogenic regeneration and have previously been shown to contribute towards the treatment of fracture nonunion [6]. Moreover, transplantation of a combination of pre-differentiated MSCs with the scaffold could improve the level of engraftment at the lesion side. Here, we demonstrated that MSCs cultured within 0.2% GrowDex® in SARSTEDT 24-well TC inserts can be driven along both adipogenic and osteogenic differentiation lineages. The xeno-free, tuneable viscosity and biocompatible nature of GrowDex® in combination with 24 well TC inserts provides a promising tool in studying adipogenic and osteogenic differentiation of MSCs in 3D and could provide valuable insights in developing new strategies towards regeneration of bone injuries, non-union or potentially critical size defects.

References

1. Pittenger, M. F., et al. (1999). „Multilineage potential of adult human mesenchymal stem cells.“ *Science* 284(5411): 143-147.
2. Dominici, M., et al. (2006). „Minimal criteria for defining multipotent mesenchymal stromal cells. The International Society for Cellular Therapy position statement.“ *Cytotherapy* 8(4): 315-317.
3. Frith, J. E., et al. (2010). „Dynamic three-dimensional culture methods enhance mesenchymal stem cell properties and increase therapeutic potential.“ *Tissue Eng Part C Methods* 16(4): 735-749.
4. Lo, Y.-P., et al. (2016). „Three-dimensional spherical spatial boundary conditions differentially regulate osteogenic differentiation of mesenchymal stromal cells.“ *Scientific Reports* 6: 21253.
5. Jung, J. P., et al. (2016). „Heterogeneous Differentiation of Human Mesenchymal Stem Cells in 3D Extracellular Matrix Composites.“ *BioResearch Open Access* 5(1): 37-48.

6. Bajada, S., et al. (2007). „Successful treatment of refractory tibial non-union using calcium sulphate and bone marrow stromal cell implantation.“ The Journal of Bone and Joint Surgery. British volume 89-B(10): 1382-1386.
7. Reznikov, N., et al. (2016). „A materials science vision of extracellular matrix mineralization.“ Nature Reviews Materials 1: 16041.
8. Azoidis, I., J., et al. (2017). „Three-dimensional cell culture of human mesenchymal stem cells in nanofibrillar cellulose hydrogels.“ MRS Communications: 1-8.

For additional product or technical information, please e-mail us at marketing@sarstedt.com or visit www.sarstedt.com.

SARSTEDT AG & Co. KG

Sarstedtstraße 1
D-51588 Nümbrecht
Germany

www.sarstedt.com
info@sarstedt.com

Worldwide Support Offices

 Australia Tel: +61 8 8349 6555 Fax: +61 8 8349 6882 info.au@sarstedt.com	 Austria Tel: +43 2236 616 82 Fax: +43 2236 620 93 info.at@sarstedt.com	 Belgium Tel: +32 3 541 7692 Fax: +32 3 541 8103 info.be@sarstedt.com	 Brazil Tel: +55 11 4152 2233 Fax: +55 11 4152 3198 info.br@sarstedt.com	 Canada Tel: +1 514 328 6614 Toll free: 1 888 727 7833 Fax: +1 514 328 9391 info.ca@sarstedt.com	 China Tel: +86 21 5062 0181 Fax: +86 21 5058 0700 info.cn@sarstedt.com
 Croatia Tel: +385 95 36 77 030 Fax: +385 14 96 10 75 info.hr@sarstedt.com	 Czech Republic Tel: +420 281 021 491 Fax: +420 281 021 495 info.cz@sarstedt.com	 Finland Tel: +358 9 374 1044 Fax: +358 9 374 1176 info.fi@sarstedt.com	 France Tel: +33 3 84 31 95 95 Fax: +33 3 84 31 95 99 info.fr@sarstedt.com	 Germany Telefon +49 22 93 305-0 Telefax +49 22 93 305-3450 Service (0800) 0 83 305-0 info@sarstedt.com	 Greece Tel: +30 210 6038 274 Fax: +30 210 6038 276 info.gr@sarstedt.com
 Hungary Tel: +36 1 383 1216 Fax: +36 1 383 1213 info.hu@sarstedt.com	 Ireland Tel: +353 53 91 44922 Fax: +353 53 91 44998 info.ie@sarstedt.com	 Italy Tel: +39 045 8510 114 Fax: +39 045 8510 118 info.it@sarstedt.com	 Japan Tel: +81 3 5215 5400 Fax: +81 3 5215 6400 info.jp@sarstedt.com	 Latvia Tel: +371 6 731 0386 Fax: +371 6 704 0723 info.lv@sarstedt.com	 Mexico Tel: +52 55 8501 1577 Fax: +52 55 8501 1578 info.mx@sarstedt.com
 Netherlands Tel: +31 76 501 7550 Fax: +31 76 501 7626 info.nl@sarstedt.com	 Norway Tel: +47 64 856 820 Fax: +47 64 856 821 info.no@sarstedt.com	 Poland Tel: +48 22 722 0543 Fax: +48 22 722 0795 info.pl@sarstedt.com	 Portugal Tel: +351 21 915 6010 Fax: +351 21 915 6019 info.pt@sarstedt.com	 Russia Tel: +7 495 937 5228 Fax: +7 495 937 5228 info.ru@sarstedt.com	 Slovakia Tel: +421 2 682 45 933 Fax: +421 2 682 45 934 info.sk@sarstedt.com
 Spain Tel: +34 93 846 4103 Fax: +34 93 846 3978 info.es@sarstedt.com	 Sweden Tel: +46 42 19 84 50 Fax: +46 42 19 84 59 info.se@sarstedt.com	 Switzerland Tel: +41 81 750 1880 Fax: +41 81 750 1899 info.ch@sarstedt.com	 Turkey Tel: +90 216 290 18 65 Fax: +90 216 290 18 64 info.tr@sarstedt.com	 United Arab Emirates Tel: +971 4 3888 080 Fax: +971 4 3888 282 info.ae@sarstedt.com	 United Kingdom Tel: +44 116 2359 023 Fax: +44 116 2366 099 info@sarstedt.co.uk
 USA Tel: +1 800 257 5101 Tel: +1 828 465 4000 Fax: +1 828 465 4003 customerservice@sarstedt.us					