

AST-2 tube stripper




Automatic and rapid stripping
of blood bag tubes

- Tube stripper with high-performance motor
- Mains or battery-operated
- Versions with and without auto-stop function
- Table-top version available with auto-stop function and foot switch

AST-2 tube stripper

AST-2 Basic model

For easy and fatigue-free stripping of blood bag tubes




AST-2S with auto-stop function

As basic model but with automatic motor stop at the blood bag



AST-2F table-top device

As basic model but with auto-stop function, bench-top plate and foot pedal operation



Technical data

Model	AST-2	AST-2S	AST-2F
Order No.	87.7709.029	87.7709.025	87.7709.035
Right or left-handed	✓	indicate when ordering	-
Overload protection	✓	✓	✓
Forward and reverse	✓	✓	✓
Auto-stop function	-	✓	✓
Foot pedal	-	-	✓
Mains operation	✓	✓	✓
Battery-powered operation	✓	✓	✓
Weight	690 g	690 g	1.480 g
High-performance motor	9 V 20 W	9 V 20 W	9 V 20 W
Dimensions of device with mains adapter (LxWxH in mm)	225x39x62	225x47x62	250x140x71
Transport case	✓	✓	-
Accessories required	AST-2, AST-2S, AST-2F		
Universal mains adapter INPUT: 100 V–240 V AC +/- 10%, 50–60 Hz, 700 mA OUTPUT: 7.5 V DC 3.3 A	87.7709.033		
Battery pack NiMh 7.2 V 700 mA	87.7709.031		
Universal charger for battery pack PRI: 100 V–240 V AC, 50–60 Hz 17 VA SEC: 4.8 V–12 V DC max. 800 mA 9.6 VA	87.7709.032		

SARSTEDT, Inc.
1025 St. James Church Road
P.O. Box 468
Newton, NC 28658 - 0468
Toll free: +1 800 257 5101
Fax: +1 828 465 4003
customerservice.us@sarstedt.com
www.sarstedt.com