

# Packaging change 50 ml tubes



## Benefits:

- Convenient pressure cap for hygienic opening and closing
- Tamper-evident closure for secure sealing
- Tubes packed in straight rows for convenient removal

## Ordering information

Art. N°.	Name	Cap colour	Certification	Pcs./inner pack	Pcs./outer box
62.547.254	Tubes, 50 ml, 114 x 28 mm, conical, PP	■	Free from DNA, DNase/RNase, pyrogens and endotoxins, non-cytotoxic, sterile	25	300
62.547.274	Tubes, 50 ml, 114 x 28 mm, conical, PP	☒	Free from DNA, DNase/RNase, pyrogens and endotoxins, non-cytotoxic, sterile	25	300
62.548.004	Tubes, 50 ml, 114 x 28 mm, conical, PP	■	Free from DNA, DNase/RNase, pyrogens and endotoxins, non-cytotoxic	25	300

Available from 06/2021