



S-Monovette® Needle / S-Monovette® Safety-Needle

Blood Collection Needle for S-Monovette® System

S-Monovette® Needle

Image	Order number	Ø/Length	Description	Color code	Packaging
	85.1160	0.9 mm/38 mm	20G x 1½"		500/case • 50/inner box individual sterile
	85.1162	0.8 mm/38 mm	21G x 1½"		500/case • 50/inner box individual sterile
	85.1440	0.7 mm/38 mm	22G x 1½"		500/case • 50/inner box individual sterile
	85.1372	0.9 mm/25 mm	20G x 1"		500/case • 50/inner box individual sterile
	85.1373	0.8 mm/25 mm	21G x 1"		500/case • 50/inner box individual sterile
	85.1441	0.7 mm/25 mm	22G x 1"		500/case • 50/inner box individual sterile

S-Monovette® Safety-Needle

Image	Order number	Ø/Length	Description	Color code	Packaging
	85.1160.200	0.9 mm/38 mm	20G x 1½"		500/case • 50/inner box individual sterile
	85.1162.200	0.8 mm/38 mm	21G x 1½"		500/case • 50/inner box individual sterile
	85.1162.600	0.8 mm/38 mm	21G x 1½" metal poor		500/case • 50/inner box individual sterile
	85.1440.200	0.7 mm/38 mm	22G x 1½"		500/case • 50/inner box individual sterile
	85.1372.200	0.9 mm/25 mm	20G x 1"		500/case • 50/inner box individual sterile
	85.1373.200	0.8 mm/25 mm	21G x 1"		500/case • 50/inner box individual sterile
	85.1441.200	0.7 mm/25 mm	22G x 1"		500/case • 50/inner box individual sterile












Multifly® Needle / Safety-Multifly® Needle

Blood Collection Needle for S-Monovette® System

Multifly® Needle

Image	Order number	Ø/Length	Description	Color code	Tubing length	Packaging
	85.1637.035	0.9 mm/ 19 mm	20G x ¾"		200 mm	480/case • 120/inner box individual sterile
	85.1638.035	0.8 mm/ 19 mm	21G x ¾"		200 mm	480/case • 120/inner box individual sterile
	85.1640.035	0.6 mm/ 19 mm	23G x ¾"		200 mm	480/case • 120/inner box individual sterile
	85.1637.005	0.9 mm/ 19 mm	20G x ¾"		80 mm	480/case • 120/inner box individual sterile
	85.1638.005	0.8 mm/ 19 mm	21G x ¾"		80 mm	480/case • 120/inner box individual sterile
	85.1640.005	0.6 mm/ 19 mm	23G x ¾"		80 mm	480/case • 120/inner box individual sterile
	85.1642.005	0.5 mm/ 19 mm	25G x ¾"		80 mm	480/case • 120/inner box individual sterile
	85.1638	0.8 mm/ 19 mm	21G x ¾" detachable guide sleeve		200 mm	480/case • 120/inner box individual sterile
	85.1638.001	0.8 mm/ 19 mm	21G x ¾" without multi-adaptor		240 mm	480/case • 120/inner box individual sterile
						

Safety-Multifly® Needle

Image	Order number	Ø/Length	Description	Color code	Tubing length	Packaging
	85.1637.235	0.9 mm/ 19 mm	20G x ¾"		200 mm	480/case • 120/inner box individual sterile
	85.1638.235	0.8 mm/ 19 mm	21G x ¾"		200 mm	480/case • 120/inner box individual sterile
	85.1640.235	0.6 mm/ 19 mm	23G x ¾"		200 mm	480/case • 120/inner box individual sterile
	85.1642.235	0.5 mm/ 19 mm	25G x ¾"		200 mm	480/case • 120/inner box individual sterile
	85.1637.205	0.9 mm/ 19 mm	20G x ¾"		80 mm	480/case • 120/inner box individual sterile
	85.1638.205	0.8 mm/ 19 mm	21G x ¾"		80 mm	480/case • 120/inner box individual sterile
	85.1640.205	0.6 mm/ 19 mm	23G x ¾"		80 mm	480/case • 120/inner box individual sterile
	85.1642.205	0.5 mm/ 19 mm	25G x ¾"		80 mm	480/case • 120/inner box individual sterile
	85.1638.200	0.8 mm/ 19 mm	21G x ¾" detachable guide sleeve		200 mm	480/case • 120/inner box individual sterile
	85.1640.200	0.6 mm/ 19 mm	23G x ¾" detachable guide sleeve		200 mm	480/case • 120/inner box individual sterile
