

Go for Gold!

The S-Monovette® –
the gold standard
in blood collection




Venous blood is the most important sample material for answering medical questions. And the basis for performing diagnostics smoothly and precisely is optimal specimen quality. The S-Monovette® from SARSTEDT has set standards in this field since its launch:

No other blood collection system provides such optimal specimen quality, regardless of the vein conditions – reliable blood collection for a reliable diagnosis.

Set standards in blood collection

The S-Monovette®





Saves time – Save time when it matters.



Safe – Protection where it is needed.



First-rate patient care – The well-being of your patients is non-negotiable for us.




Sustainable – An investment in a better future.



Everything from a single source – Achieve more together. Let's work together!



Ready to use immediately thanks to the integrated needle holder (multi-adapter).



The slim, ergonomic design of the needle holder ensures a flat puncturing angle.

Special needle guard minimises the risk of needlestick injuries.



SARSTEDT

SARSTEDT AG & Co. KG

P.O. Box 12 20

D-51582 Nümbrecht

Phone: +49 2293 305 0

Fax: +49 2293 305 3992

export@sarstedt.com

www.sarstedt.com



Further information, videos
and trial devices:
s-monovette.sarstedt.com

If you have any questions:

We'd be happy to help!

Visit our website: www.sarstedt.com



SARSTEDT