

# Cell culture products

by SARSTEDT



**READY.  
SET.  
GROW!**

# SARSTEDT: Excellence in cell culture



Try it out for free and  
with no obligation!

cellculture.  
[sarstedt.com/en](https://www.sarstedt.com/en)

Set new standards in cell cultivation with SARSTEDT! Our innovative product line, consisting of flasks, dishes, plates and a wide range of complements, provides many different growth surfaces for all in vitro conditions and is tailored to the requirements of most cell types. Our precise production and highest quality standards pave the way for excellent science.

The SARSTEDT colour codes make it easy to find the right receptacle at a glance. Put your trust in the scientific expertise and product portfolio of SARSTEDT.

## Cell culture dishes

- SARSTEDT SUREGrip – Outer ring enables a **secure grip through haptic feedback**
- Clear arrows on the base and lid for **uniform positioning**
- **Gas exchange** through air vents in the lid
- **Safe stacking** thanks to distinct rings
- Dishes in  $\varnothing$  35 mm and  $\varnothing$  60 mm with a **grid for cloning experiments**

### The SARSTEDT colour code

indicates the three surface properties of the culture vessels that create optimal conditions for cell growth.

**Red** – The standard surface for adherent cells

**Yellow** – The SARSTEDT Cell+ surface for difficult adherent cells and for serum-reduced/serum-free conditions of cultivated cells

**Green** – The SARSTEDT suspension culture surface

Find the right surface for your cell culture here:



## Cell culture flasks

- **Accessible** for pipettes and scrapers
- Large **labelling fields**
- **Engraved and printed** graduated scales
- **Good stability** with stacking rim
- **Reduction of contamination risk** due to angled neck and anti-drip rim

## Cell culture plates

- Alphanumeric labelling on the edge and between the wells for **better orientation** when filling
- **Reduced risk of contamination** when pipetting thanks to free-standing wells
- **Secure hold** due to non-slip side grids
- Transparent side walls allow **visual inspection** of the medium
- Integrated air vents and condensation rings in the lid for **gas exchange and minimisation of evaporation**



We support the entire  
cell culture workflow with an  
extensive product portfolio.

For an overview, see [cellculture.sarstedt.com/en](http://cellculture.sarstedt.com/en)  
**Let's grow! SARSTEDT**

## **SARSTEDT AG & Co. KG**

Sarstedtstraße 1  
D-51588 Nümbrecht

Phone: +49 2293 305 0  
Fax: +49 2293 305 3992

[export@sarstedt.com](mailto:export@sarstedt.com)  
[www.sarstedt.com](http://www.sarstedt.com)