

Cell culture products

by SARSTEDT



**READY.
SET.
GROW!**



Life Science

SARSTEDT: Excellence in cell culture



Try it out for free and
with no obligation!

[shop.sarstedt.us/
cellculture](https://shop.sarstedt.us/cellculture)

Set new standards in cell cultivation with SARSTEDT! Our innovative product line, consisting of flasks, dishes, plates and a wide range of complementary products, provides many different growth surfaces for all in vitro conditions and is tailored to the requirements of most cell types. Our precise production and highest quality standards pave the way for excellent results.

SARSTEDT color codes make it easy to find the right vessel at a glance. Put your trust in the scientific expertise and product portfolio of SARSTEDT.

Cell culture dishes

- SARSTEDT SUREGrip – Outer ring enables a **secure grip at a touch**
- Clear arrows on the base and lid for **uniform positioning**
- **Gas exchange** through air vents in the lid
- **Safe stacking** with distinct rings
- Dishes in \varnothing 35 mm and \varnothing 60 mm with **grid for cloning experiments**

Cell culture flasks

- **Accessible** for pipettes and scrapers
- Large **labelling fields**
- **Molded and printed** graduations
- **Good stability** with stacking rim
- **Reduction of contamination risk** due to canted neck and anti-drip rim

Cell culture plates

- Alphanumeric labeling on the edge and between the wells for **better orientation** when filling
- **Reduced risk of contamination** when pipetting via free-standing wells
- **Secure hold** due to non-slip side grids
- Transparent side walls allow **visual inspection** of the medium
- Integrated air vents and condensation rings in the lid for **gas exchange and to minimize evaporation**

The SARSTEDT colour code

indicates three different cell culture surfaces designed to create optimal conditions for cell growth.

Red – The standard surface for adherent cells

Yellow – The SARSTEDT Cell+ surface for difficult adherent cells and for serum-reduced/serum-free conditions of cultivated cells

Green – The SARSTEDT suspension culture surface

Find the right surface for your cell culture here:



We support the entire
cell culture workflow with an
extensive product portfolio.

For an overview, see shop.sarstedt.us/cellculture
Let's grow! SARSTEDT

SARSTEDT, Inc.

1025 St. James Church Road
P.O. Box 468
Newton, NC 28658 - 0468

Toll free: +1 800 257 5101

customerservice.us@sarstedt.com
www.sarstedt.com