



### **Product description**

Order number 87.7709.031  
Product description Battery Battery for AST-2 tube stripper

### **Product characteristics**

Application Obtaining blood and blood products  
Supply voltage 7.2 V  
Battery power 700 mAh

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 1  
Type of smallest subpackaging case  
Piece(s) / inner box 1  
Piece (s) / outer case 1  
Piece(s) / pallet 432  
Depth of box 41 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification Battery

Page 2

Width of box	41 mm
Height of box	102 mm
Depth of case	188 mm
Width of case	188 mm
Height of case	100 mm
Case volume	0.0035 cbm
Weight of product	0.109 kg
Weight of case	0.26 kg
EAN of inner box	4038917141178
EAN case	4038917279697

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**