



Product description

Order number	74.5330.520
Product description	UD 3plus, Urine drainage system, 2 l, with holder, sterile, tube length: 120 cm, 1 piece(s)/blister
PZN number	16086831

Product characteristics

Colour of print/label	blue
Design	with holder
Application	Urine drainage
Volume of work	2 l

Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	10.1 mm
Length of tube	120 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polyvinylchloride (PVC)
Colour of product transparent

Purity & certification

Product category Medical device
CE certified CE supplier
Purity standard sterile
Sterilisation Gas sterilisation
Batched yes

Packaging

Minimum order qty. 20
Type of smallest blister
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 20
Piece(s) / pallet 840
Depth of case 395 mm
Width of case 295 mm
Height of case 220 mm
Case volume 0.0256 cbm
Weight of product 0.179 kg
Weight of case 3.58 kg
EAN of inner box 6955824024340
EAN case 6955824024197

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

