



### Product description

Order number	74.5485.005
Product description	Urine drainage system, 2 l, with anti-reflux valve, with T drainage valve, sterile, tube length: 100 cm, 1 piece(s)/blister
PZN number	16356199

### Product characteristics

Colour of print/label	blue
Design	with T drainage valve
Application	Urine drainage
Volume of work	2 l
Anti-reflux valve	yes

### Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	10.1 mm
Length of tube	100 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine drainage system, 2 l, with anti-reflux valve,  
with T drainage valve, sterile**

**Page 2**

**Material & colours**

Product material Polyvinylchloride (PVC)  
Colour of product transparent

**Purity & certification**

Product category Medical device  
CE certified CE supplier  
Purity standard sterile  
Sterilisation Gas sterilisation  
Batched yes

**Packaging**

Minimum order qty. 50  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 1  
Piece (s) / outer case 50  
Piece(s) / pallet 2400  
Depth of case 395 mm  
Width of case 295 mm  
Height of case 290 mm  
Case volume 0.0338 cbm  
Weight of product 0.0844 kg  
Weight of case 4.62 kg  
EAN of inner box 6955824024418  
EAN case 6955824024036

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**