




Product Specification
Urine measuring system, 2 l, with anti-reflux valve, exchange bag 2,000 ml, sterile



Product description

Order number	74.5681.200
Product description	Urine measuring system, 2 l, with anti-reflux valve, exchange bag 2,000 ml, sterile, tube length: 120 cm, 1 piece(s)/blister

Product characteristics

Colour of print/label	blue
Design	exchange bag 2,000 ml
Application	Urine drainage
Volume of work	500 ml
Anti-reflux valve	yes

Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	10.1 mm
Length of tube	120 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



Product Specification
Urine measuring system, 2 l, with anti-reflux valve, exchange bag 2,000 ml, sterile

Page 2

Material & colours

Product material Polyvinylchloride (PVC)
Colour of product transparent

Purity & certification

Product category Medical device
CE certified CE supplier
Purity standard sterile
Sterilisation Gas sterilisation
Batched yes

Packaging

Minimum order qty. 10
Type of smallest blister
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 10
Piece(s) / pallet 280
Depth of case 580 mm
Width of case 290 mm
Height of case 235 mm
Case volume 0.0395 cbm
Weight of product 0.3364 kg
Weight of case 4.04 kg
EAN of inner box 6955824024425
EAN case 6955824024043

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT