



**Product Specification**  
**Urine measuring system exchange bag, 3 l, with drainage valve, sterile**



**Product description**

Order number	74.5685.002
Product description	Urine measuring system exchange bag, 3 l, with drainage valve, sterile, 1 piece(s)/blister

**Product characteristics**

Colour of print/label	blue
Design	with drainage valve
Application	Urine drainage
Volume of work	3 l

**Size**

Sample volume	3 l
Max. volume	3 l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



### **Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

### **Purity & certification**

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

### **Packaging**

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	4050
Depth of case	395 mm
Width of case	245 mm
Height of case	200 mm
Case volume	0.0194 cbm
Weight of product	0.0473 kg
Weight of case	3.22 kg
EAN of inner box	6955824024463
EAN case	6955824024081

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

