



### **Product description**

Order number	74.5767.003
Product description	UroStretch, size: L, material: PL
Aid number	15.25.11.1034
PZN number	8406960

### **Material & colours**

Product material	polyester (PL)
Colour of product	white

### **Purity & certification**

Product category	Medical device
CE certified	CE supplier
Batched	yes

### **Packaging**

Minimum order qty.	10
Type of smallest subpackaging	case
Piece(s) / inner box	10
Piece (s) / outer case	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**UroStretch, size: L**

Piece(s) / pallet	2340
Depth of case	283 mm
Width of case	173 mm
Height of case	135 mm
Case volume	0.0066 cbm
Weight of product	0.043 kg
Weight of case	0.55 kg
EAN of inner box	5703936100121
EAN case	5703936100169

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

