




Product Specification

Urine bag, 1.5 l, with anti-reflux valve, with drainage valve

Page 1



Product description

Order number	74.5218.101
Product description	Urine bag, 1.5 l, with anti-reflux valve, with drainage valve, tube length: 90 cm, 150 piece(s)/case
PZN number	16333689

Product characteristics

Colour of print/label	blue
Design	with drainage valve
Application	Urine drainage
Volume of work	1,5 l
Anti-reflux valve	yes

Size

Sample volume	1.5 l
Max. volume	1.5 l
Diameter of tube	6.6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Urine bag, 1.5 l, with anti-reflux valve, with drainage valve

Page 2

Material & colours

Product material Polyvinylchloride (PVC)
Colour of product transparent

Purity & certification

Product category Medical device
CE certified CE supplier
Batched yes

Packaging

Minimum order qty. 150
Type of smallest case
subpackaging
Piece(s) / inner box 150
Piece (s) / outer case 150
Piece(s) / pallet 9600
Depth of case 580 mm
Width of case 390 mm
Height of case 150 mm
Case volume 0.0339 cbm
Weight of product 0.0352 kg
Weight of case 5.8 kg
EAN of inner box 6971243660969
EAN case 6971243660693

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT