



Product description

Order number	74.5283.101
Product description	Urine bag, 1.5 l, with anti-reflux valve, tube length: 90 cm, 200 piece(s)/case
PZN number	16333778

Product characteristics

Colour of print/label	blue
Application	Urine drainage
Volume of work	1,5 l
Anti-reflux valve	yes

Size

Sample volume	1.5 l
Max. volume	1.5 l
Diameter of tube	6.6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polyvinylchloride (PVC)
Colour of product transparent

Purity & certification

Product category Medical device
CE certified CE supplier
Batched yes

Packaging

Minimum order qty. 200
Type of smallest case
subpackaging
Piece(s) / inner box 200
Piece (s) / outer case 200
Piece(s) / pallet 9600
Depth of case 580 mm
Width of case 390 mm
Height of case 150 mm
Case volume 0.0339 cbm
Weight of product 0.0325 kg
Weight of case 7 kg
EAN of inner box 6971243660952
EAN case 6971243660686

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

