



Product description

Order number	97.8710.800
Product description	SAHARA-III MAXITHERM 230 V, for thawing and warming blood products in plastic bags before transfusion, max. capacity: 6 storage bags, suitable for plasma, whole blood, erythrocyte concentrates and cryopreserved preparations

Product characteristics

Application	Warming
Supply voltage	230 V
Power consumption during operation	655 W
Frequency	50 Hz - 60 Hz

Size

Width of product	32 cm
Height of product	32.5 cm
Length of product	49.3 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	Medical device, CE
CE certified	CE
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Case volume	0.1713 cbm
Weight of product	14.44 kg
Weight of case	16.4 kg
EAN of inner box	4038917968614
EAN case	4038917968621

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

