



Product description

| | |
|---------------------|--|
| Order number | 97.8710.602 |
| Product description | SAHARA-TSC 115 V, for thawing and warming blood products in plastic bags before transfusion, suitable for cryopreserved stem cells, incl. adaptation compress and warming tray |

Product characteristics

| | |
|------------------------------------|---------------|
| Application | Warming |
| Supply voltage | 115 V |
| Power consumption during operation | 655 W |
| Frequency | 50 Hz - 60 Hz |

Size

| | |
|-------------------|---------|
| Width of product | 32 cm |
| Height of product | 32.5 cm |
| Length of product | 49.3 cm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|--------------------|
| Product category | Medical device, CE |
| CE certified | CE |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 1 |
| Piece(s) / pallet | 1 |
| Case volume | 0.1639 cbm |
| Weight of product | 21 kg |
| Weight of case | 23.75 kg |
| EAN of inner box | 4038917968577 |
| EAN case | 4038917968584 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

