



Product description

Order number	97.8710.620
Product description	Warming shell, Replacement part for storing stem cell preparations

Product characteristics

Application	Warming
-------------	---------

Size

Width of product	18 cm
Length of product	27 cm

Material & colours

Product material	aluminium
------------------	-----------

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Case volume	0.0096 cbm
Weight of product	2.45 kg
Weight of case	2.95 kg
EAN of inner box	4028339720525
EAN case	4028339720518

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

