



Product description

Order number	19.930.125
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 125 µl, (LxØ): 110 x 2.05 mm, 250 piece(s)/container

Product characteristics

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	70 IU/ml blood
Application	blood gas analysis

Size

Sample volume	125 µl
Diameter	2.05 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polyethylene therephthalate (PET)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Latex-free yes
Batched yes

Packaging

Minimum order qty. 250
Type of smallest container
subpackaging
Piece(s) / inner box 250
Piece (s) / outer case 16000
Piece(s) / pallet 1280000
Depth of case 378 mm
Width of case 288 mm
Height of case 180 mm
Case volume 0.0196 cbm
Weight of product 0.0003 kg
Weight of case 6.64 kg
EAN of inner box 4038917234283
EAN case 4038917234276

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

