



Product description

Order number	19.930.175
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 175 µl, (LxØ): 100 x 2.3 mm, 200 piece(s)/container

Product characteristics

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	70 IU/ml blood
Application	blood gas analysis

Size

Sample volume	175 µl
Diameter	2.3 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polyethylene terephthalate (PET)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Latex-free yes
Batched yes

Packaging

Minimum order qty. 200
Type of smallest container
subpackaging
Piece(s) / inner box 200
Piece (s) / outer case 12800
Piece(s) / pallet 921600
Depth of case 333 mm
Width of case 303 mm
Height of case 190 mm
Case volume 0.0192 cbm
Weight of product 0.0003 kg
Weight of case 5.5 kg
EAN of inner box 4038917234344
EAN case 4038917234337

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

