



Product Specification
Quick-release cap 1.60?2.10 mm for blood gas capillaries 100 µl / 125 µl, white



Product description

| | |
|---------------------|---|
| Order number | 65.935.205 |
| Product description | Quick-release cap 1.60?2.10 mm for blood gas capillaries 100 µl / 125 µl, white, blood gas analysis, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------|---------------------------------------|
| Type of collection | capillary |
| Application | blood gas analysis |
| Suitable for | blood gas capillaries 100 µl / 125 µl |

Material & colours

| | |
|-------------------|--------------------------------|
| Product material | Thermoplastic elastomers (TPE) |
| Colour of product | white |

Purity & certification

| | |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Batched | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Quick-release cap 1.60?2.10 mm for blood gas capillaries 100 µl / 125 µl, white

Page 2

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 32000 |
| Piece(s) / pallet | 2560000 |
| Depth of case | 378 mm |
| Width of case | 288 mm |
| Height of case | 180 mm |
| Case volume | 0.0196 cbm |
| Weight of product | 0.0002 kg |
| Weight of case | 8.42 kg |
| EAN of inner box | 4038917234559 |
| EAN case | 4038917234542 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT