



### **Product description**

Order number 90.1090  
Product description ESR rack, suitable for S-Sedivette® , 10 measuring stations

### **Product characteristics**

Suitable for S-Sedivette®  
Number of positions 10

### **Material & colours**

Product material Polypropylene (PP)

### **Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	320
Depth of case	598 mm
Width of case	398 mm
Height of case	350 mm
Case volume	0.0833 cbm
Weight of product	0.6018 kg
Weight of case	11.75 kg
EAN of inner box	4038917085304
EAN case	4038917004008

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

