



Product description

Order number	90.1092.730
Product description	Barcode scanner, incl. USB connection, for Sediplus® S 2000 NX

Product characteristics

Suitable for	Sediplus® S 2000 NX
--------------	---------------------

Size

Width of product	85 mm
Height of product	60 mm
Length of product	150 mm

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	324
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.366 kg
Weight of case	0.47 kg
EAN of inner box	4038917927130
EAN case	4038917927147

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

