



### **Product description**

Order number	85.1015
Product description	Safety lancet, Mini, Ø needle: 28 G, light blue, penetration depth: 1.6 mm, for capillary blood collection, for puncturing finger, ear, heel, 200 piece(s)/case

### **Product characteristics**

Application	capillary blood collection
Example of use	Young children, adults: diabetic patients, oncology, determining blood sugar
Puncture site	finger, ear, heel
Penetration depth	1.6 mm
Safety product	yes

### **Size**

Blood volume	0,1 µl - 10 µl
Diameter of needle	28 G
Width of product	21 mm
Length of product	52 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Safety lancet, Mini, Ø needle: 28 G, penetration depth: 1.6 mm**

**Material & colours**

Colour of product light blue

**Purity & certification**

Product category Medical device  
CE certified CE 0124 (DEKRA)  
Purity standard sterile  
Sterilisation Gamma irradiation  
Batched yes

**Packaging**

Minimum order qty. 1200  
Type of smallest case  
subpackaging  
Piece(s) / inner box 200  
Piece (s) / outer case 1200  
Piece(s) / pallet 108000  
Depth of box 118 mm  
Width of box 118 mm  
Height of box 151 mm  
Depth of case 366 mm  
Width of case 248 mm  
Height of case 173 mm  
Case volume 0.0157 cbm  
Weight of product 0.0028 kg  
Weight of case 4.08 kg  
EAN of inner box 4038917085021  
EAN case 4038917057622

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

