



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 17.2111.050   |
| Product description | Minivette® POCT neutral Z, nominal volume: 50 µl, without preparation, white, 200 piece(s)/bag, 200 piece(s)/case |

### **Product characteristics**

|                     |                            |
|---------------------|----------------------------|
| Type of collection  | capillary                  |
| Type of preparation | without preparation        |
| Application         | capillary blood collection |
| Coding colour       | white                      |

### **Size**

|               |       |
|---------------|-------|
| Sample volume | 50 µl |
|---------------|-------|

### **Material & colours**

|                   |                        |
|-------------------|------------------------|
| Product material  | Cellulose acetate (CA) |
| Colour of product | transparent            |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Latex-free       | yes                           |
| Batched          | yes                           |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 2000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 200           |
| Piece (s) / outer case        | 2000          |
| Piece(s) / pallet             | 144000        |
| Depth of box                  | 140 mm        |
| Width of box                  | 99 mm         |
| Height of box                 | 117 mm        |
| Depth of case                 | 513 mm        |
| Width of case                 | 293 mm        |
| Height of case                | 140 mm        |
| Case volume                   | 0.021 cbm     |
| Weight of product             | 0.001 kg      |
| Weight of case                | 2.66 kg       |
| EAN of inner box              | 4038917310314 |
| EAN case                      | 4038917310307 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

