



**Product Specification**  
**Microvette® APT 500 EDTA K2E, 500 µl, cap violet, cap, round base**



**Product description**

Order number	20.2330.100
Product description	Microvette® APT 500 EDTA K2E, capillary blood collection, preparation: K2 EDTA, nominal volume: 500 µl, Ø: 13 mm, with plastic label, cap violet, cap, with pierceable membrane, round base, inner receptacle: cylindrical, non-sterile, 100 piece(s)/case

**Product characteristics**

Type of collection	capillary
Type of preparation	K2 EDTA
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	transparent/violet
Cap	membrane cap
Closure type	cap
Application	capillary blood collection
Base shape	round base
Volume of work	250 µl - 500 µl
Shape of inner tube	cylindrical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Microvette® APT 500 EDTA K2E, 500 µl, cap violet, cap, round base**

**Page 2**

**Size**

Sample volume	500 µl
Diameter	13 mm
Length excluding cap	75 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polyethylene (PE)
Colour of cap	violet

**Purity & certification**

Satisfies the requirement	P 650, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	non-sterile
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	54000
Depth of box	166 mm
Width of box	176 mm
Height of box	68 mm
Depth of case	358 mm
Width of case	197 mm
Height of case	186 mm
Case volume	0.0131 cbm
Weight of product	0.006 kg
Weight of case	3.17 kg
EAN of inner box	4038917946032
EAN case	4038917946049

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**