




Product Specification
Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent, sterile



Product description

Order number	80.623
Product description	Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon, assembled, cap: brown, (LxØ) without cap: 101 x 16.5 mm, outer diameter of cap: 19.5 mm, sterile, 500 piece(s)/bag

Product characteristics

Cap	Screw cap with stool spoon
Closure type	screw cap
Design	spoon
Base shape	round base

Size

Outer diameter of cap	19.5 mm
Diameter	16.5 mm
Length of product	107 mm
Length excluding cap	101 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent, sterile

Page 2

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Stem material	Polystyrene (PS)
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	22000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0.0306 cbm
Weight of product	0.0068 kg
Weight of case	3.89 kg
EAN of inner box	4038917082983
EAN case	4038917001960

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT