



**Product Specification**  
**Faeces tube, with spoon, screw cap, (LxØ): 76 x 20 mm, transparent, sterile**



**Product description**

Order number	80.734.001
Product description	Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon, assembled, cap: brown, (LxØ) without cap: 76 x 20 mm, sterile, 100 piece(s)/bag

**Product characteristics**

Cap	Screw cap with stool spoon
Closure type	screw cap
Design	spoon
Base shape	flat-conical base

**Size**

Diameter	20 mm
Length of product	77 mm
Length excluding cap	76 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Faeces tube, with spoon, screw cap, (LxØ): 76 x 20 mm, transparent, sterile**

**Page 2**

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Stem material	Polystyrene (PS)
Tube material	Polypropylene (PP)

### Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	20000
Depth of case	598 mm
Width of case	398 mm
Height of case	165 mm
Case volume	0.0393 cbm
Weight of product	0.0069 kg
Weight of case	4.06 kg
EAN of inner box	4038917119801
EAN case	4038917016810

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**