



### Product description

Order number	80.9924.027
Product description	Faeces tube, with blade, white, screw cap, cap: brown, (LxØ) without cap: 55 x 44 mm, with paper label, 250 piece(s)/bag

### Product characteristics

Label/ Print	with paper label
Colour of print/label	black
Closure type	screw cap
Design	blade

### Size

Diameter	44 mm
Length of product	57 mm
Length excluding cap	55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Faeces tube, with blade, screw cap, (LxØ): 55 x 44 mm, white**

**Page 2**

**Material & colours**

Colour of product	white
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	8000
Depth of case	458 mm
Width of case	408 mm
Height of case	373 mm
Case volume	0.0697 cbm
Weight of product	0.0137 kg
Weight of case	7.8 kg
EAN of inner box	4038917154819
EAN case	4038917039635

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-15

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**