



## Product description

Order number	93.852.169
Product description	Block Rack D17, for 50 collection tubes, Ø opening: 17 mm, format: 5 x 10, violet, material: recycled PP

## Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

## Size

Diameter of opening	17 mm
Width of product	109 mm
Height of product	45 mm
Length of product	209 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material recycled Polypropylene (recycled PP)  
Colour of product violet

### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 20  
Piece(s) / pallet 760  
Depth of case 568 mm  
Width of case 433 mm  
Height of case 120 mm  
Case volume 0.0295 cbm  
Weight of product 0.1958 kg  
Weight of case 4.62 kg  
EAN of inner box 4038917122047  
EAN case 4038917044899

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

