



**Product Specification**  
**Block Rack D13, Ø opening: 13 mm, 5 x 10,  
orange**

**Page 1**



### Product description

Order number	93.853.135
Product description	Block Rack D13, for 50 collection tubes, Ø opening: 13 mm, format: 5 x 10, orange, material: PP

### Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

### Size

Diameter of opening	13 mm
Width of product	96 mm
Height of product	45 mm
Length of product	180 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Block Rack D13, Ø opening: 13 mm, 5 x 10,  
orange**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product orange

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 16  
Piece(s) / pallet 1440  
Depth of case 398 mm  
Width of case 368 mm  
Height of case 120 mm  
Case volume 0.0176 cbm  
Weight of product 0.1687 kg  
Weight of case 3.16 kg  
EAN of inner box 4038917122566  
EAN case 4038917074063

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**