



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 93.853.200  |
| Product description | Double block rack D13, for 100 collection tubes, Ø opening: 13 mm, format: 10 x 10, grey, material: recycled PP |

### **Product characteristics**

|                        |         |
|------------------------|---------|
| Number of wells        | 100     |
| Storage capacity       | 100     |
| Format                 | 10 x 10 |
| Autoclavable (max. °C) | 121 °C  |

### **Size**

|                     |        |
|---------------------|--------|
| Diameter of opening | 13 mm  |
| Width of product    | 192 mm |
| Height of product   | 45 mm  |
| Length of product   | 180 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Double block rack D13, Ø opening: 13 mm, 10 x 10, grey**

**Page 2**

**Material & colours**

Product material recycled Polypropylene (recycled PP)  
Colour of product grey

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 8  
Piece(s) / pallet 720  
Depth of case 398 mm  
Width of case 368 mm  
Height of case 120 mm  
Case volume 0.0176 cbm  
Weight of product 0.3509 kg  
Weight of case 3.27 kg  
EAN of inner box 4038917122573  
EAN case 4038917074070

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**