



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, flat tube bottom, PP, without  
cap, 100 piece(s)/bag**

**Page 1**



### Product description

|                     |  |
|---------------------|--|
| Order number        | 60.611.311   |
| Product description | Screw cap tube, working volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, flat tube bottom, brown, material: PP, without cap, 100 piece(s)/bag, 1,000 piece(s)/carton |

### Product characteristics

|                     |  |
|---------------------|--|
| Light protection    | yes                                    |
| Base shape          | conical false bottom, flat tube bottom |
| Volume of work      | 5 ml                                   |
| Shape of inner tube | conical                                |

### Size

|                      |         |
|----------------------|---------|
| Diameter             | 15.3 mm |
| Height of product    | 92 mm   |
| Length of product    | 92 mm   |
| Length excluding cap | 92 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, flat tube bottom, PP, without  
cap, 100 piece(s)/bag**

**Page 2**

### Material & colours

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | brown              |
| Tube material     | Polypropylene (PP) |

### Purity & certification

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 30000         |
| Case volume                   | 0.046 cbm     |
| Weight of product             | 0.004 kg      |
| Weight of case                | 4.55 kg       |
| EAN of inner box              | 4038917367943 |
| EAN case                      | 4038917367936 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**