



**Product Specification**  
**Archiving cap, light blue, suitable for**  
**S-Monovette® , tubes Ø 13-16 mm**

**Page 1**



**Product description**

Order number 65.647.020  
Product description Archiving cap, with filter, light blue, suitable for S-Monovette® , tubes Ø 13-16 mm, 1,000 piece(s)/bag

**Product characteristics**

Suitable for S-Monovette® , tubes Ø 13-16 mm  
Filter yes

**Material & colours**

Product material Polypropylene (PP)  
Colour of product light blue  
Colour of cap light blue

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Archiving cap, light blue, suitable for**  
**S-Monovette® , tubes Ø 13-16 mm**

**Page 2**

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.001 kg
Weight of case	5.37 kg
EAN of inner box	4038917250085
EAN case	4038917250078

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**