



Product description

| | |
|---------------------|---|
| Order number | 75.570.001 |
| Product description | Container, max. volume: 100 ml, (LxØ): 49 x 70 mm, Ø opening: 70 mm, transparent, cap: transparent, material: PS, snap-on lid, without print, non-sterile, 50 piece(s)/tube bag |

Product characteristics

| | |
|------------------|---------------|
| Label/ Print | without print |
| Graduation | no |
| Closure type | snap-on lid |
| Volume of work | 100 ml |
| Shape of opening | round |

Size

| | |
|---------------------|--------|
| Sample volume | 100 ml |
| Max. volume | 100 ml |
| Diameter | 70 mm |
| Diameter of opening | 70 mm |
| Height of product | 49 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Container, max. volume: 100 ml, (LxØ): 49 x 70 mm, PS, transparent, without print

Page 2

Length of product 49 mm
Length excluding cap 49 mm

Material & colours

Product material Polystyrene (PS)
Colour of product transparent
Colour of cap transparent

Purity & certification

Product category no medical device | IVD
Purity standard non-sterile
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging tube bag
Piece(s) / inner box 50
Piece (s) / outer case 1000
Piece(s) / pallet 45000
Depth of case 375 mm
Width of case 295 mm
Height of case 315 mm
Case volume 0.0348 cbm
Weight of product 0.0027 kg
Weight of case 3.1 kg
EAN of inner box 4038917974370
EAN case 4038917974387

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT