



### **Product description**

Order number	76.571.001
Product description	Snap-on lid, for 100-ml container 75.570.001, material: PET, transparent, 50 piece(s)/tube bag

### **Product characteristics**

Closure type	snap-on lid
Suitable for	100-ml container 75.570.001

### **Material & colours**

Product material	Polyethylene therephthalate (PET)
Colour of product	transparent

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	tube bag
Piece(s) / inner box	50
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Depth of case	375 mm
Width of case	295 mm
Height of case	315 mm
Case volume	0.0348 cbm
Weight of product	0.001 kg
Weight of case	1.4 kg
EAN of inner box	4038917974462
EAN case	4038917974479

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

