



Product Specification
**S-Monovette® , 25 ml, for Formalin system, (LxØ):
97 x 25 mm, with paper label**

Page 1



Product description

Order number 09.1705.001
Product description S-Monovette®, up to 25 ml, for Formalin system, (LxØ): 97 x 25 mm, with paper label, 50 piece(s)/bag

Product characteristics

Type of preparation without preparation
Label/ Print with paper label
Colour of print/label black
Closure type screw cap
Application Tissue fixation
Suitable for Formalin system
Volume of work 25 ml

Size

Max. volume 25 ml
Diameter 25 mm
Length of product 97 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® , 25 ml, for Formalin system, (LxØ):
97 x 25 mm, with paper label

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	12000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0141 kg
Weight of case	4.01 kg
EAN of inner box	4038917335287
EAN case	4038917335270

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT