



### **Product description**

Order number	73.1052.001
Product description	Insert, suitable for tubes and S-Monovette® Ø 15 mm, material: PS, length: 32 mm, transparent, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Suitable for	tubes and S-Monovette® Ø 15 mm
Base shape	skirted conical base

### **Size**

Diameter of opening	11 mm
Length of product	32 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      2000  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      500  
Piece (s) / outer case                      2000  
Piece(s) / pallet                            156000  
Depth of case                              428 mm  
Width of case                                323 mm  
Height of case                              130 mm  
Case volume                                0.018 cbm  
Weight of product                        0.0008 kg  
Weight of case                              1.98 kg  
EAN of inner box                         4038917124874  
EAN case                                      4038917019507

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

