



Product description

Order number	73.1055
Product description	RIA tube, suitable for gamma counters, material: PP, 0.5 ml, (LxØ): 22 x 8 mm, transparent, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	gamma counters
Base shape	conical base

Size

Max. volume	0.5 ml
Diameter	8 mm
Diameter of opening	29 mm
Height of product	22 mm
Length of product	22 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
RIA tube, suitable for gamma counters,
transparent

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 10000
Type of smallest bag
subpackaging
Piece(s) / inner box 1000
Piece (s) / outer case 10000
Piece(s) / pallet 480000
Depth of case 388 mm
Width of case 378 mm
Height of case 215 mm
Case volume 0.0315 cbm
Weight of product 0.0005 kg
Weight of case 5.05 kg
EAN of inner box 4038917085175
EAN case 4038917003865

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT