



### **Product description**

Order number	73.646
Product description	Sample tube, suitable for auto-analyzer, material: PS, 3.5 ml, (LxØ): 38 x 14 mm, transparent, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Suitable for	auto-analyzer
Base shape	conical base

### **Size**

Max. volume	3.5 ml
Diameter	14 mm
Diameter of opening	8 mm
Height of product	38 mm
Length of product	38 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, suitable for auto-analyzer,  
transparent**

**Page 2**

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 5000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 1000  
Piece (s) / outer case 5000  
Piece(s) / pallet 120000  
Depth of case 398 mm  
Width of case 398 mm  
Height of case 415 mm  
Case volume 0.0657 cbm  
Weight of product 0.0018 kg  
Weight of case 9.62 kg  
EAN of inner box 4038917083096  
EAN case 4038917002066

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**