



### **Product description**

Order number 73.650.500  
Product description Sample tube, suitable for Coulter Counter, cap enclosed, material: PS, (LxØ): 60 x 34 mm, 250 piece(s)/bag

### **Product characteristics**

Graduation no  
Closure type push cap  
Suitable for Coulter Counter  
Base shape flat base

### **Size**

Diameter 34 mm  
Diameter of opening 12 mm  
Length of product 60 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Closure material	Low Density Polyethylene (LD-PE)

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	18000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0077 kg
Weight of case	4.51 kg
EAN of inner box	4038917157841
EAN case	4038917157834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

