



Product description

Order number	73.662.500
Product description	Scintillation vials, material: HD-PE, 20 ml, (LxØ): 58 x 27 mm, 500 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Base shape	flat base

Size

Max. volume	20 ml
Diameter	27 mm
Diameter of opening	17 mm
Height of product	58 mm
Length of product	58 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material High Density Polyethylene (HD-PE)
Closure material High Density Polyethylene (HD-PE)

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 500
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 500
Piece(s) / pallet 21000
Depth of case 380 mm
Width of case 370 mm
Height of case 235 mm
Case volume 0.033 cbm
Weight of product 0.006 kg
Weight of case 4 kg
EAN of inner box 4038917115728
EAN case 4038917014861

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

