



Product description

Order number	73.663
Product description	Sample tube, suitable for Boehringer ES 300® and ES 600®, material: PS, 2 ml, (LxØ): 40 x 12 mm, transparent, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Boehringer ES 300® and ES 600®
Base shape	conical base

Size

Max. volume	2 ml
Diameter	12 mm
Diameter of opening	10 mm
Height of product	40 mm
Length of product	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Sample tube, suitable for Boehringer ES 300® and ES 600®, transparent

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 5000
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 5000
Piece(s) / pallet 180000
Depth of case 398 mm
Width of case 398 mm
Height of case 300 mm
Case volume 0.0475 cbm
Weight of product 0.0012 kg
Weight of case 6.72 kg
EAN of inner box 4038917083195
EAN case 4038917002158

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-15

This document was prepared by EDP support and is valid without signature.



SARSTEDT